

sonicare

Sonic electric toothbrush - Trial

DiamondClean

5 modes

2 brush heads

Glass charger

USB travel charger



HX9342/03



Ultimate cleaning, superior results

Superior plaque removal plus powerful whitening

Only Philips Sonicare DiamondClean provides the performance your patients need and the features they want. DiamondClean removes up to 7x more plaque than a manual brush and eliminates surface stains to whiten smiles in just one week.

Provides a superior clean

- DiamondClean toothbrush head for Sonicare's best whitening
- Patented Sonicare toothbrush technology

Guides you to follow dental recommendations

• 5 modes: Clean, Whitening, Polish, Sensitive, Gum Care

Designed around you

- $\bullet \ \ \mathsf{Dual} \ \mathsf{charging} \ \mathsf{system}$
- · Illuminated display
- Three-week battery life

Highlights

DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond shaped bristles to effectively yet gently remove surface stains and scrape away plaque. Removes stain build up from coffee, tea, red wine, cigarettes and tobacco in just 1 week

Five brushing modes



Includes the stain-removing White mode, clinically proven to whiten teeth in as little as one week.

Sonic technology



Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

Charging system



Innovative glass charger and station for home use, plus portable charging case with laptopready USB and wall plug for travel.

Illuminated display



When on, unique backlight indicates selected mode. When off, display becomes invisible for a clean look.

Rechargeable battery



When fully charged, DiamondClean has enough battery life for three weeks of regular

Specifications

Modes

Clean: For exceptional everyday clean Gum Care: Gently massages gums Polish: Brighten and polish your teeth Sensitive: Gentle teeth and gum cleaning White: Removes surface stains

Items included

Handles: 1 DiamondClean

Brush heads: 1 DiamondClean standard, 1

DiamondClean compact Glass charger: 1

Travel case: USB travel charger

Design and finishing

Color: Ceramic white

Cleaning performance

Speed: Up to 62000 brush movement/min Health benefits: Improves gum health in only two weeks

Performance: Removes up to 100% more

olaque

Whitening benefits: Whitens teeth up to 2

shades

Timer: SmarTimer and Quadpacer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Light shows battery status

Handle: Slim ergonomic design Display: Illuminated display

Technical specifications

Battery: Rechargeable Battery type: Lithium ION

Operating time (full to empty): 3 week

battery life

Service

Warranty: 2-year limited warranty

Power

Voltage: 110-220 V

© 2019 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2019-07-24 Version: 17.1.1

EAN: 00 07502 00206 11

www.philips.com



^{* *}based on two periods of two-minute brushings per day, on clean mode