

PHILIPS

sonicare

Rechargeable powered
interdental cleaner

AirFloss Pro

2 nozzles

250 ml mouthwash



HX8432/12



Healthier gums in 2 weeks, guaranteed*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro is the easiest way to effectively clean between teeth. AirFloss Pro can be used with mouthwash or water and is clinically proven to be as effective as floss for gum health.**

Improves oral health

- Up to 99.9% plaque removal***
- Clinically proven to be as effective as floss for gum health**
- Helps prevent cavities between teeth

Superior technology

- Triple-burst customisable settings
- High-performance nozzle

The easiest way to effectively clean between teeth

- An easy way to start a healthy habit

Achieve a 30-second clean

- Point. Press. Clean!

Easy to use

- Technology that cleans with a touch
- Guidance tip ensures correct placement

Highlights

Better gum health




Philips Sonicare AirFloss Pro/Ultra nozzle is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.***

High-performance nozzle




The new AirFloss Pro/Ultra high-performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Air and micro-droplet technology



AirFloss utilises unique air and micro-droplet technology to send a burst of pressurised air and water or mouthwash between your teeth. The one-touch activation button is as simple as a mouse click.

Helps prevent cavities




By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

Guidance tip



The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

Triple-burst customisation




AirFloss Pro nozzle's powerful spray can be customised to your liking, with single, double or triple bursts for each press of the activation button.

Easy-fill reservoir



Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

Specifications

Items included

AirFloss Pro handle: 1
Nozzle: 2
Charger: 1
Mouthwash: 237 ml BreathRx

Design and finishing

Colour: Pink

Cleaning performance

Clinically proven results: As effective as floss, for gum-health improvement, for reducing gum bleeding, for plaque reduction

Performance: Removes up to 99.9% plaque**

Health benefits: Cleans entire mouth in 60 secs, Improves gum health in 2 weeks

Customisable settings: Auto-burst or manual burst, Single, double, triple burst

For best results: Change nozzle every 6 months

Ease of use

Handle: Slim ergonomic design

Battery indicator: Handle LEDs show charge level

Battery Life: Up to 2 weeks

Nozzle system: Easily snaps on and off

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Charge time: 24 hours to full charge

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty



* or your money back

** when used in conjunction with a manual toothbrush and anti-microbial mouthwash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under the Support tab for further details.

*** In a lab study; actual in-mouth results may vary.