



HX8372/01

Healthier gums in 2 weeks, guaranteed*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro is the easiest way to effectively clean between teeth. AirFloss Pro can be used with mouthwash or water and is clinically proven to be as effective as floss for gum health.**

Improves oral health

- · Clinically proven to be as effective as floss for gum health**
- Up to 99.9% plaque removal***
- · Helps prevent cavities between teeth

Superior technology

- Air and micro-droplet technology
- High-performance nozzle
- · Triple-burst customisable settings

The easiest way to effectively clean between teeth

- Takes just 60 seconds: Point, Press, Clean!
- · Total daily confidence with a fresh and clean mouth
- · An easy way to start a healthy habit

Highlights

Better gum health



Philips Sonicare AirFloss Pro is clinically proven to improve gum health as much as floss.** Helps to improve gum health in as little as two weeks.

Up to 99.9% plaque removal




Removes up to 99.9% of plaque from treated areas.***

Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro helps prevent cavities from forming in the spaces between your teeth.

Micro-droplet technology



Our clinically proven results are possible thanks to our unique technology, which combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

High-performance nozzle



The new AirFloss Pro high-performance nozzle amplifies the power of our air and microdroplet technology to make it more effective and efficient than ever.

Triple-burst customisation



AirFloss Pro's powerful spray can be customised to your liking, with single, double or triple bursts for each press of the activation button.

Cleans in just 60 seconds



With the AirFloss Pro, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple) and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Fresh and clean with mouthwash



Fill the reservoir on the handle with either mouthwash or water, and then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

Specifications

Items included

AirFloss Pro handle: 1

Nozzle: 2 Charger: 1

Cleaning performance

Clinically proven results: As effective as floss, for gum-health improvement, for reducing gum bleeding, for plaque reduction

Performance: Removes up to 99.9% of plaque Health benefits: Cleans entire mouth in 60 secs, Improves gum health in 2 weeks Customisable settings: Auto-burst or manual burst, Single, double, triple burst For best results: Change nozzle every 6 months

Ease of use

Handle: Slim ergonomic design

Battery indicator: Handle LEDs show charge

level

Battery Life: Up to 2 weeks

Nozzle system: Easily snaps on and off

Technical specifications

Battery: Rechargeable Battery type: Lithium ION

Charge time: 24 hours to full charge

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty

Design and finishing

Colour: White with grey accents

© 2019 Koninklijke Philips N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

Issue date 2019-07-23 Version: 4.0.1

EAN: 08 71010 37150 09

www.philips.com



- * * when used in conjunction with a manual toothbrush and anti-microbial mouthwash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under the Support tab for further details.
- * * * In a lab study; actual in-mouth results may vary.
- * or your money back