

# PHILIPS

## sonicare

Rechargeable sonic  
toothbrush

Elite

RS 900 series



HX6942/10



## Redefine clean

with FlexCare features

Innovative technology delivers better results. Now there's a guaranteed way to better oral health. It's the new Philips Sonicare FlexCare that adapts to your oral care needs.

### ProResults brush head

- Designed for a better fit with your mouth

### FlexCare features

- 3 flexible brushing modes, 2 personalised cleaning routines
- Go Care: 1-minute brushing cycle for a quick clean
- Max care: 3 minutes for an extended care routine

### UV sanitiser

- All-in-one sanitising, charging and storing solution

# Highlights

## ProResults brush head



ProResults brush head's contour fit, wider bristle range and broader sweeping motion provide better tooth coverage.

## FlexCare features



3 flexible brushing modes, 2 personalised cleaning routines.

## Go Care

One-minute brushing cycle in the Clean mode, for a quick clean.

## Max care

3 minutes for an extended care routine.

## UV sanitiser



All-in-one sanitising, charging and storing solution.

## Features

### Rechargeable battery

**Smartimer:** 2-min. timer helps to ensure recommended brushing

**Quadpacer:** Interval timer

**Easy-Start:** Gently increases power over first 14 uses

**Deluxe recharge gauge:** Light signal helps track battery charge

## Technical specifications

**Voltage:** Multi-Voltage Charger V

**Speed:** 31,000 brush strokes per minute

**Battery type:** Lithium Ion rechargeable battery

**Operating time (full to empty):** Fully charged should provide 20 2-min. brushings

## Weight and dimensions

**F-box dimensions:** 6.75 x 9.25 x 3.75

(WxHxD) in

**F-box weight:** 1.93 lb

**A-box dimensions:** 11.825 x 7.25 x 10

(LxWxH) in

**A-box weight:** 6.3 lb

## Design and finishing

**Chrome front panel**

**Colour(s):** Medical Blue

**Material handle:** ABS/PA, TPE

**Material brush head:** PCTA, TPE, with nylon bristles

## Items Included

**UV sanitiser:** 1

**Deluxe soft travel case:** 1

**FlexCare handle:** 1

**Standard ProResults brush head:** 1

**Compact ProResults brush head:** 1

**Additional brush head colour code rings:**  
3 pcs

