

# PHILIPS

## sonicare

Compact sonic  
toothbrush heads

InterCare

1-pack

Standard size

Click-on

Cleaning in between teeth

HX9011/08



## Reaches deeper between teeth

Designed to improve gum health in two weeks

Philips Sonicare InterCare brush heads are designed for patients who want a deeper clean between teeth. Clinically proven to remove up to 7x more plaque between teeth than a manual toothbrush after four weeks of use.\*

### Provides a superior clean

- Designed to reach deeper and go further
- Engineered to maximize sonic motion

### Designed around you

- Reminder bristles
- Works with any Philips Sonicare click-on toothbrush

### Proven to improve oral health

- Part of a better oral healthcare routine



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Sonic technology




Philips Sonicare brush heads are vital to our core technology of high-frequency, high-amplitude brush movements performing over 31,000 brush strokes per minute. This dynamic cleaning action drives fluid deep between teeth and along the gumline, for a thorough yet gentle cleaning.

## Reach deeper and go further



This Philips Sonicare toothbrush head features extra-long, high-density bristles to target hidden plaque caught deep between teeth and other hard-to-clean areas.

## Reminder bristles



At first glance it may not be obvious, but brush heads lose stiffness and gradually wear down over months of normal use. Our blue reminder

bristles fade to white and help patients recognize when it's time for a replacement. For optimal results, replace brush heads every three months.

## Click-on design



Philips Sonicare brush heads click on and off the brush handle for a secure fit and easy maintenance and cleaning. Fits all Philips Sonicare toothbrush handles except: PowerUp and Essence.

## Better oral healthcare



Like all authentic Philips Sonicare branded brush heads, this brush head is safe on teeth and gums. Each brush head has been quality tested for exceptional performance and durability.

# Specifications

## Items included

**Brush heads:** 1 InterCare compact

## Compatibility

**Brush head system:** Click-on

**Suitable for these models:** 2 Series plaque control, 3 Series gum health, DiamondClean, EasyClean, FlexCare, FlexCare Platinum, FlexCare+, for Kids, HealthyWhite, HealthyWhite+

## Design and finishing

**Color:** White

**Bristle stiffness feel:** Regular

**Reminder bristles:** Blue bristle color fade away

**Size:** Compact

## Health benefits

**Gum health:** Improves gum health in 2 weeks

**Plaque removal:** Removes up to 7x more plaque\*

## Quality and performance

**Replacement:** Every 3 months

**Tested:** for optimal usage



[asimpleswitch.com](https://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

