

PHILIPS

sonicare

Sonic electric
toothbrush

ProtectiveClean
4300

Built-in pressure sensor

1 cleaning mode

1 x BrushSync feature



HX6888/63



Cleaner teeth. Gently does it.

Remove up to 7 x more plaque than a manual

Feel the difference of a gentle clean with our Pressure Sensor while removing up to 7 x more plaque than a manual toothbrush

Innovative technology

- Improves gingival health in just 2 weeks
- Connects smart brush handle and smart brush heads
- A gentle nudge to go nice and easy
- Philips Sonicare's advanced sonic technology
- Patients get the most from their brush heads

Designed around your patient


- A great way for patients to get used to their toothbrush
- Encouragement to brush thoroughly

Improves patients' oral health

- Fillings, crowns, or sensitive areas? Not to worry.

Highlights

Improves gingival health



The i InterCare brush head features extra-long, densely packed, high-quality bristles to target deep interdental plaque. It is clinically proven to reduce gingival bleeding and inflammation in just two weeks.

Safe and gentle




You can be sure of a safe brushing experience: our sonic technology is suitable for use with braces, fillings, crowns and veneers, and helps prevent cavities and improve gum health.

Philips Sonicare technology



Philips Sonicare's advanced sonic technology pulses water between your teeth, and its brush strokes break plaque up and sweep it away for an exceptional daily clean.

Smart replacement reminders



All brush heads wear out over time. But we've got that covered too: our smart brush-recognition technology tracks how long patients have been using their brush head for and how hard they've been brushing. When it's time to replace it, a light on the handle lets them know. That way, they can be sure that their brush head is consistently doing a stellar job.

Too much pressure?




Patients may not notice if they are brushing too hard, but their toothbrush will. If patients start taking it out on their gums, the toothbrush will make a pulsing sound as a reminder to ease off on the pressure.

BrushSync technology



A microchip-enabled technology that detects and synchronises the smart brush head with the smart handle. The smart handle and smart brush head pair is a powerful combination that enables Smart replacement reminders.

Easy-start programme



Starting something new can take getting used to. That's why our easy-start programme gives patients the option of a gradual, gentle increase in brushing power the first 14 times they use their new toothbrush.

Helpful timers



Our QuadPacer lets patients know when they've spent just the right amount of time cleaning each part of their mouth, while our Smartimer tells them when they've brushed for the recommended two minutes.

Specifications

Modes

Clean: For exceptional everyday cleaning

2 intensities: Low, High

Items included

Handle: 1 ProtectiveClean

Brush heads: 1 i InterCare standard

Charger: 1

Design and finishing

Colour: Mid-Blue

Cleaning performance

Speed: Up to 62,000 brush movements/min

Performance: Removes up to 7x more plaque*

Health benefits: Helps reduce cavities

Pressure feedback: Vibrates handle to alert user

Timer: SmartTimer and Quadpacer

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

BrushSync Replacement Reminder: Always know when to, replace brush heads

Ease of use

Handle compatibility: Easy click-on brush heads

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Brushing time: Up to 2 weeks

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 2 weeks

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

