

# PHILIPS

## sonicare

Sonic electric  
toothbrush

EasyClean

1 mode

1 brush head

HX6511/43



## Better plaque removal

Sonicare's unique dynamic cleaning action gently and effectively reaches deep between teeth and along the gumline.

### **Provides a superior clean**

- Removes up to 7 x more plaque than a manual toothbrush
- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth

### **Proven to improve oral health**

- Philips Sonicare toothbrush helps whiten teeth
- Clinically proven safe and gentle

### **Offers a customised brushing experience**

- Easy-start programme to acclimate to the Sonicare experience

### **Guides you to follow dental recommendations**

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Helps whiten teeth



Philips Sonicare's best teeth whitening: this toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

## Removes up to 7 x more plaque



Removes up to 7 times more plaque from hard-to-reach places than a manual toothbrush.

## Sonic technology



The Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

## Angled brush head neck



The uniquely angled neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

## Safe and gentle




Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

## Smartimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

## Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

## Easy-start programme



Gently increases power over first 14 uses to ease you into the Philips Sonicare electric toothbrush experience



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Modes

**Clean:** For exceptional everyday cleaning

## Items included

**Handles:** 1 EasyClean

**Brush heads:** 1 DiamondClean standard

**Charger:** 1

## Design and finishing

**Colour:** White

## Cleaning performance

**Speed:** Up to 31,000 brush strokes/min

**Performance:** Removes up to 7x more plaque\*

**Health benefits:** Helps improve gum health

**Whitening benefits:** Whitens teeth by up to 2 shades

**Timer:** Smartimer and Quadpacer

## Ease of use

**Brush head system:** Easy click-on brush heads

**Brushing time:** Up to 2 weeks

**Battery indicator:** Light shows battery status

**Handle:** Slim ergonomic design

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** Up to 2 weeks

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

