

PHILIPS

sonicare

Sonic electric
toothbrush

HealthyWhite

3 modes

2 brush heads



HX6732/02



Whiter, healthier teeth

The ultimate toothbrush for whiter teeth

Everybody loves a bright white smile. Let the natural whiteness of your teeth shine through. Sonicare HealthyWhite electric toothbrush HX6732/02 is proven to remove everyday stains in just two weeks, with the use of the Clean & White Mode.



Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Removes up to 2x more plaque than a manual toothbrush
- Dynamic cleaning action drives fluid between teeth

Proven to improve oral health

- Naturally whiter teeth with patented sonic technology
- Philips Sonicare toothbrush helps whiten teeth
- Sensitive Mode: Gentle teeth and gum cleaning

Offers a customised brushing experience

- Easy-start programme to acclimate to the Sonicare experience
- Clean and White Mode: proven to remove stains
- Clean Mode: Standard 2-minute mode


Highlights

Patented sonic technology



The Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gum line.

Naturally whiter teeth




The dynamic cleaning action of this toothbrush and the direct broader contact with each tooth, has been proven to remove everyday stains, to help keep teeth naturally whiter.

Smartimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Quadpacer



30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on, resulting in more consistent cleaning throughout your mouth

Easy-start programme



Gently increases power over first 14 uses to ease you into the Philips Sonicare electric toothbrush experience

Clean and White Mode




2 minutes of Clean mode with an additional 30 seconds of White mode to focus on your visible front teeth. Removes everyday stains, such as coffee, tea, tobacco and red wine. Whitens teeth by 2 shades in just 2 weeks.*

Clean Mode


For exceptional everyday cleaning. Helps maintain whitening treatments.

Helps whiten teeth




Philips Sonicare's best teeth whitening: this toothbrush helps to remove and reduce stains on your teeth for a brighter smile.

Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.

Gentle teeth and gum cleaning



Philips Sonicare toothbrush with Sensitive Mode: Gentle, yet thorough cleaning for sensitive gums

Patented Sonic technology



The Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

Specifications

Modes

Clean: For exceptional everyday cleaning

Sensitive: Gentle teeth and gum cleaning

Clean and White: 2 min Clean + 30 sec White

Items included

Handles: 1 HealthyWhite

Brush heads: 1 ProResults compact, 1 ProResults standard

Charger: 1

Travel case: 1

Design and finishing

Colour: Frost White

Cleaning performance

Speed: Up to 31,000 brush strokes/min

Performance: Removes up to 90% of stains*

Health benefits: Helps improve gum health

Whitening benefits: Whitens teeth by up to 2 shades

Timer: Smartimer and Quadpacer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Light shows battery status

Handle: Slim ergonomic design, Rubber grip for easy handling

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 3 weeks**

Power

Voltage: 110–220 V

Service

Warranty: 2-year limited warranty

