

PHILIPS

& TensRelief 4 Self-adhesive electrodes

PulseRelief

Fits on PR3840, PR3093,
PR3094

4 electrodes

Patented hydrogel technology

Easy magnet connection

PR3820



For an active life with less pain

Recognized pain relief therapy by professionals

The Philips self-adhesive electrodes have patented hydrogel technology to minimize skin irritation. There are 4 in a pack and recommended for up to 20 uses. The electrodes are easy to snap on the TENS Units with a magnet connection.

Patented hydrogel technology

- Patented hydrogel technology to minimize skin irritation

Optimised for multiple use

- Recommended to use for multiple applications up to 20 times

Magnet for easy connection

- Easy magnet connection to connect TENS Units with electrodes

Highlights

Minimize skin irritation

Our electrodes use a high quality patented hydrogel technology to minimize the effects of skin irritation. The electrodes are flexible so they can be applied on irregular body contours.

Recommended use up to 20 times

To be used for multiple applications on different body parts. The number of times the electrodes can be used is dependent on skin

type and usage. The electrodes can easily be stored in their original plastic package for next time use.

Easy snap-on connection

The TENS units are easily connected to the self-adhesive electrodes by a magnet. First, apply the self-adhesive electrodes to your body and then just snap-on the TENS units via a magnet connection.

Specifications

What's in the box

Self-adhesive electrodes: 4 pcs

Self-adhesive electrodes size: 50 x 50 mm

Multistick® hydrogel

Plastic seal bag

Composition of hydrogel

Glycerin, non-tallow: (Cas. No.: 56-81-5)

Water: (Cas. No.: 28062-44-4)

Poly(acrylate) Co-polymer: (Cas. No.: 7732-18-5)

Country of origin

USA

