

PHILIPS

sonicare

Rechargeable sonic
toothbrush - Trial

FlexCare Platinum

3 modes, 3 intensities

1 brush head

With pressure sensor

UV Brush Head Sanitizer



HX9141/23



Deep cleaning

Removes up to 7x more plaque between teeth*

With nine unique brushing options, an intuitive pressure sensor and advanced brush head technology, Philips Sonicare FlexCare Platinum goes beyond traditional cleaning to deliver exceptional plaque removal and improved gum health.

Offers a customized brushing experience

- Customize your mode and intensity

Provides an exceptional clean for teeth and gums

- Clinically proven to improve gum health in two weeks
- Removes up to 7x more plaque between teeth
- Designed to reach deeper and go further
- Patented sonic technology for better oral health

Fits your lifestyle

- Three week battery life between charges
- Handy travel case and charger

Helps you meet recommended brushing guidelines

- Alerts you when you're brushing too hard
- Timers to help encourage thorough brushing



asimpleswitch.com

Highlights

9 total brushing experiences




FlexCare Platinum lets you choose from three brushing modes – Clean, White or Deep Clean – and three intensity settings – Low, Medium or High – to create a brushing experience that’s perfectly tailored to you. Clean mode for an exceptional everyday clean (2 minutes), White mode for stain removal (2:30 minutes) and Deep Clean mode for an invigorating deep clean (3 minutes).

Improves gum health in two wks




With a special Gum Care mode for focused cleaning to help reduce gingival bleeding and inflammation, FlexCare Platinum is clinically proven to improve gum health in only two weeks.

Reaches deep between teeth



FlexCare Platinum reaches deep between teeth to provide superior plaque removal. It is clinically proven to remove 7x more plaque between teeth than a manual toothbrush.

Intuitive Pressure Sensor



FlexCare Platinum features a unique, intuitive Pressure Sensor to help reduce aggressive brushing. The handle gently vibrates to let you know when you’re applying too much pressure.

Reach deeper and go further



Philips Sonicare InterCare toothbrush head features extra-long, high-density bristles to target hidden plaque caught deep between teeth and other hard-to-clean areas.

Long battery life



When fully charged, your FlexCare Platinum has enough battery life for three weeks of regular use.

Travel case and charger



It’s easy to take your FlexCare Platinum wherever you need it with an easy-to-pack travel case and compact travel charger.

Patented sonic technology



Unlike a traditional electric or power toothbrush, Sonicare features a patented sonic technology. Sonicare’s unique dynamic cleaning action drives fluid between teeth, gently and effectively removing plaque interdentially and along the gum line.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Modes

Clean: For exceptional everyday clean

White: Removes surface stains

Deep Clean: For an invigorating deep clean

3 intensities: High, Low, Medium

Items included

Handles: 1 FlexCare Platinum

Brush heads: 1 InterCare standard

Charger: 1

UV sanitizer

Travel case: 1

Design and finishing

Color: Gray with brushed white panel

Cleaning performance

Speed: Up to 31,000 brush strokes/min

Performance: Removes up to 7x more plaque*

Health benefits: Improves gum health in only two weeks

Whitening benefits: Helps naturally remove stains

Timer: Smartimer and Quadpacer

Pressure sensor: Vibrates handle to alert user

Ease of use

Brush head system: Easy snap-on brush heads

Brushing time: Up to 3 weeks**

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Display: Illuminated display

Technical specifications

Battery: Rechargeable

Operating time (full to empty): Up to 3 weeks

Battery type: Lithium ION

Power

Voltage: 110–220 V

Service

Warranty: 2-year limited warranty

