

# PHILIPS

## sonicare

Sonic electric  
toothbrush – Trial

FlexCare Platinum

5 modes

3 brush heads



HX9143/08



## More innovation, less plaque between teeth Effective interdental cleaning for teeth and gums

Philips Sonicare FlexCare Platinum goes beyond traditional cleaning to remove up to 7x more plaque between teeth than a manual toothbrush. An intuitive pressure sensor and nine brushing combinations allow for personalized care and guidance.

### **Provides a superior clean**

- Patented Sonicare toothbrush technology

### **Designed around you**

- InterCare brush head

### **Guides you to follow dental recommendations**

- Intuitive pressure sensor

### **Offers a customized brushing experience**

- Customize your mode and intensity

### **Proven to improve oral health**

- Quadpacer
- Smarttimer

# Highlights

## Sonic technology




Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

## InterCare brush head



Extra-long, high-density bristles reach deeper between teeth to help remove significantly more plaque\* interproximally. Ideal for patients with bleeding gums or increased risk of gum inflammation.

## Intuitive pressure sensor



If too much pressure is applied, the handle gently pulsates to tell patients they are brushing too hard to guide proper technique. In a study of U.S. patients, 7 out of 10 aggressive brushers tested reduced their brushing pressure when using the pressure sensor.

## 9 total brushing experiences




FlexCare Platinum lets you choose from three cleaning modes — Clean, White or Gum Care — and three intensity settings — Low, Medium or High — to create a brushing experience that's perfectly tailored to you.

## Smarttimer



Two-minute timer helps ensure dental professional recommended brushing time.

## Quadpacer



Thirty-second interval timer encourages thorough brushing of each quadrant of the mouth.

# Specifications

## Modes

**3 intensities:** Low, Medium, High

**Clean:** For exceptional everyday clean

**Gum Care:** Gently massages gums

**White:** Removes surface stains

## Items included

**Brush heads:** 1 DiamondClean compact, 1

InterCare compact, 1 Sensitive compact

**Charger:** 1

**Travel case:** 1

**Handles:** 1 FlexCare Platinum

## Design and finishing

**Color:** White w/ metallic front panel

## Cleaning performance

**Speed:** Up to 62000 brush movement/min

**Performance:** Removes up to 7x more plaque\*

**Health benefits:** Improves gum health in only two weeks

**Whitening benefits:** Helps naturally remove stains

**Timer:** SmartTimer and Quadpacer

**Pressure sensor:** Vibrates handle to alert user

## Ease of use

**Battery indicator:** Light shows battery status

**Brush head system:** Easy click-on brush heads

**Brushing time:** Up to 3 weeks\*\*

**Handle:** Slim ergonomic design

**Display:** Illuminated display

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Operating time (full to empty):** 3 week battery life

## Power

**Voltage:** 110-220 V

## Service

**Warranty:** 2-year limited warranty

