

PHILIPS

sonicare

Sonic electric
toothbrush

EasyClean

1 mode

2 brush heads

HX6512/02



Better plaque removal

The switch from a manual toothbrush made easy

The unique dynamic cleaning action of the Philips Sonicare electric toothbrush gently and effectively reaches deep between teeth and along the gumline.

Guides you to follow dental recommendations

- Two-minute timer helps ensure recommended brushing time
- Quadpacer interval timer encourages thorough brushing

Provides a superior clean

- Sonicare dynamic cleaning action drives fluid between teeth
- Angled brush head neck for better reach to back teeth
- Removes up to 2x more plaque than a manual toothbrush
- Dynamic cleaning action drives fluid between teeth

Proven to improve oral health

- Clinically proven safe and gentle
- Philips Sonicare toothbrush helps whiten teeth

Offers a customised brushing experience

- Easy-start program to acclimate to the Sonicare experience



asimpleswitch.com

Highlights

Safe and gentle



Your Philips Sonicare Electric Toothbrush is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

Sonic technology



Philips Sonicare electric toothbrush's unique dynamic action gently and effectively reaches deep between teeth and along the gumline.

SmarTimer




2-minute timer on this Philips Sonicare electric toothbrush helps ensure dental professional recommended brushing time

Angled brush head neck



The uniquely angled brush head neck of this toothbrush head makes it easier to reach back teeth, removing plaque in those hard-to-reach areas.

Easy-start program



Gently increases power over first 14 uses to ease into the Philips Sonicare electric toothbrush experience

Quadpacer




30 seconds interval timer indicates when you have completed each quadrant of your mouth and signals for you to move on resulting in more consistent cleaning throughout your mouth

Helps whiten teeth



Helps to remove and reduce stains on your teeth for a brighter smile.

Removes up to 2x more plaque



Patented Sonic Technology removes up to 2x more plaque than a manual toothbrush.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Modes

Clean: For exceptional everyday clean

Items included

Handles: 2 EasyClean

Brush heads: 2 ProResults standard

Charger: 1

Design and finishing

Color: Glacier Green

Cleaning performance

Speed: Up to 31,000 brush movements/min

Performance: Up to 2X more plaque removal than a manual toothbrush

Health benefits: Protects gum health

Whitening benefits: Helps whiten teeth

Timer: SmartTimer and Quadpacer

Ease of use

Brush head system: Easy click-on brush heads

Brushing time: Up to 2 weeks

Battery indicator: Light shows battery status

Handle: Slim ergonomic design

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): Up to 2 weeks

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

