PHILIPS



Sonic electric toothbrush - Trial

DiamondClean

5 modes 4 brush heads USB travel charger with Deep Clean mode





Ultimate cleaning, superior results

Superior plaque removal plus powerful whitening

Only Philips Sonicare DiamondClean provides the performance your patients need and the features they want. DiamondClean removes up to 10x more plaque than a manual brush and eliminates surface stains to whiten smiles in just one week.

Provides a superior clean

- DiamondClean toothbrush head for Sonicare's best whitening
- AdaptiveClean brush head
- New Deep Clean: Effortlessly tackle trouble spots
- Patented Sonicare toothbrush technology

Guides you to follow dental recommendations

• 5 modes: Clean, White, Sensitive, Gum Care, Deep Clean

Designed around you

- Illuminated display
- Three-week battery life
- Charging system

DiamondClean brush head



DiamondClean toothbrush heads, Philips Sonicare's best whitening brush head, have diamond shaped bristles to effectively yet gently remove surface stains and scrape away plaque. Removes stain build up from coffee, tea, red wine, cigarettes and tobacco in just 1 week.

AdaptiveClean brush head



Soft, flexible sides made of rubber allow the bristles to conform to the unique shape of patients' teeth and gums to provide up to 4x more surface contact*. The brush head's flexibility helps the bristles gently track along the gumline and absorb excessive brushing pressure. Enhanced sweeping motion delivers a unique in-mouth feel and superior cleaning along the gumline and in hard-to-reach areas*. Proven to remove 10x more plaque than a manual toothbrush in hard-to-reach areas*.

Five brushing modes



Includes the stain-removing White mode, clinically proven to whiten teeth in as little as one week.

Deep Clean mode



Our all-new Deep Clean mode delivers an exceptional clean in hard-to-reach places. The unique motion of the brush head in this mode plus an extended three-minute brushing time makes it the ideal choice when you want to pay extra attention to trouble spots.

Charging system



Portable charging case with laptop-ready USB and wall plug for travel.

Illuminated display



When on, unique backlight indicates selected mode. When off, display becomes invisible for a clean look.

Rechargeable battery



When fully charged, DiamondClean has enough battery life for three weeks of regular use.

Sonic technology



Drives fluid between the teeth and along the gumline for a powerful yet gentle clean.

Specifications

Modes

Clean: For exceptional everyday clean White: Removes surface stains Sensitive: Gentle teeth and gum cleaning Gum Care: Gently massages gums Deep Clean: For an invigorating deep clean

Items included

Handles: 1 DiamondClean Brush heads: 1 DiamondClean compact, 1 DiamondClean standard, 1 InterCare standard, 1 AdaptiveClean Travel case: USB travel charger

Design and finishing

Color: Ceramic white

Cleaning performance

Speed: Up to 62000 brush movement/min Health benefits: Improves gum health in only two weeks Performance: Removes up to 10x more plaque* Whitening benefits: Whitens teeth up to 2 shades Timer: SmarTimer and Quadpacer

Ease of use

Brush head system: Easy click-on brush heads Brushing time: Up to 3 weeks** Battery indicator: Light shows battery status Handle: Slim ergonomic design Display: Illuminated display

Technical specifications

Battery: Rechargeable Battery type: Lithium ION Operating time (full to empty): 3 week battery life

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

Specifications are subject to change without notice. Trademarks are the

property of Koninklijke Philips N.V. or

their respective owners.

Issue date 2019-07-24 Version: 4.0.1





* than a manual toothbrush

* *based on two periods of two-minute brushings per day, on clean mode