

# PHILIPS

## sonicare

Rechargeable Sonic  
Toothbrush



HX9552



## Custom Care System

### 4 Ways of Customized Cleaning in One Product

No two people are the same – especially when it comes to oral care. So why not use a power toothbrush that can adapt to your different needs and still deliver superior results.

#### 2 Brush Head Sizes and 2 Cleaning Speeds

- Standard for a complete clean, compact for precision clean
- Thorough everyday clean, massages & stimulates gums

#### Customizable Cleaning

- Quadpacer® 30-second interval timer encourages even brushing
- Easy-start® eases users into the Sonicare® experience

#### Superior Cleaning Technology

- Dynamic fluid action reaches deep between teeth
- Removes plaque in hard to reach areas
- Clinically proven safe and gentle
- Healthier Gums

# Highlights

## Safe and Gentle

Your Sonicare is safe to use on: Braces (brush heads wear out sooner when used on braces), Dental restorations (fillings, crowns, veneers) and Periodontal pockets.

## Healthier Gums

Clinically proven to reduce gingivitis and help shrink periodontal pockets.

# Specifications

## Additional Features

**Exclusive patented sonic technology:** Dynamic fluid action helps clean between teeth

**Smartimer®:** Dentists' recommend a 2 minute brushing time

**2 Contoured brush head sizes:** Standard for complete clean, Compact for precision

**2 Speeds:** Choose the speed most comfortable for you

**Easy-Start®:** Gently increases power over 14 uses

**Programmable Quadpacer®:** 30-second interval timer for an even clean

**Deluxe recharge gauge:** Green tracks battery charge, yellow for recharge

**Ergonomic tapered handle:** Non-slip grip designed to improve handling

## Technical Specifications

**Operating Time (Full to Empty):** 40

**Charging time:** 24 hour(s)

**Battery Rechargeable:** NiCd

**Voltage:** 220-240

**Power consumption:** 3 W

**Speed (Top):** Up to 31,000 movements per minute

## Accessories

**HX7001:** Standard Elite single pack brush head

**HX7002:** Standard Elite double pack brush head

**HX7011:** Mini Elite single pack brush head

**HX7012:** Mini Elite double pack brush head

## Weight and dimensions

**F-box dimensions:** 17.2 x 9.5 x 23.5 in

**F-box weight:** 680 lb

**A-box weight:** 4.08 lb

**A-box dimensions:** 35.1 x 29.7 x 24.1 in

**A-box number of pieces:** 6 pcs

