

PHILIPS

sonicare

Rechargeable powered
interdental cleaner

AirFloss Ultra

2 nozzles



HX8372/03



Healthier gums in 2 weeks, guaranteed*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Ultra is the easiest way to effectively clean between teeth. AirFloss Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.**

Improves oral health

- Up to 99.9% plaque removal***
- Clinically proven to be as effective as floss for gum health**
- Helps prevent cavities between teeth

Superior technology

- Air and micro-droplet technology
- High-performance nozzle
- Triple-burst customisable settings

The easiest way to effectively clean between teeth

- Total daily confidence with a fresh and clean mouth
- An easy way to start a healthy habit
- Takes just 60 seconds: Point, Press, Clean!


Highlights

Up to 99.9% plaque removal



Removes up to 99.9% of plaque from treated areas.**

Better gum health




Philips Sonicare AirFloss Pro nozzle is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Ultra helps prevent cavities from forming in the spaces between your teeth.

Micro-droplet technology



Our clinically proven results are possible as a result of our unique technology, which combines air and mouth rinse or water to powerfully yet gently clean between teeth and along the gumline.

High-performance nozzle




The new AirFloss Ultra high-performance nozzle amplifies the power of our air and micro-droplet technology to make it more effective and efficient than ever.

Triple-burst customisation




AirFloss Ultra's powerful spray can be customised to your liking, with single, double or triple bursts for each press of the activation button.

Fresh and clean with mouthwash



Fill the reservoir on the handle with either mouthwash or water, and then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit.

Cleans in just 60 seconds



With the AirFloss Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts, or press and release for manual burst mode.

Specifications

Items included

AirFloss Pro handle: 1
Nozzle: 2
Charger: 1

Cleaning performance

Clinically proven results: As effective as floss, for gum-health improvement, for reducing gum bleeding, for plaque reduction
Performance: Removes up to 99.9% of plaque
Health benefits: Improves gum health in 2 weeks, Cleans entire mouth in 60 secs
Customisable settings: Auto-burst or manual

burst, Single, double, triple burst
For best results: Change nozzle every 6 months

Ease of use

Handle: Slim ergonomic design
Battery indicator: Handle LEDs show charge level
Battery Life: Up to 2 weeks
Nozzle system: Easily snaps on and off

Design and finishing

Colour: Black

Technical specifications

Battery: Rechargeable
Battery type: Lithium ION
Charge time: 24 hours to full charge

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty

