PHILIPS



Sonicare AirFloss

AirFloss

Rechargeable



HX8181

An easier way to floss

Gently and effectively improves oral health

If flossing isn't a regular part of your oral care routine, you're not reaching the bacteria that can build up between teeth causing infection. Sonicare AirFloss helps you to get a deep clean everyday by taking the hassle out of flossing.

Achieve a 30-second clean

• Simple one-button use; as easy as a mouse-click

Proven improvement of oral health

• Up to 99.9% plaque removal

Simple: point & press

- Easy-to-fill reservoir
- Guidance tip for easy placement
- Microburst technology cleans at the touch of a button
- Slim, angled nozzle for easy to even the hard-to-reach areas
- Slim, easy to hold handle
- Fits standard Sonicare chargers
- Charging after two weeks

HX8181/02

Highlights

Easy-to-fill reservoir



Easy-to-fill reservoir uses only a tiny amount of water for less mess. Fill with mouthwash for enhanced oral health benefits.

Fits Sonicare chargers



You only need one charger to charge your Sonicare products.

Guidance tip



Simply glide the guidance tip along the gumline until you feel it settle between the teeth. Provides easy positioning even on those hardto-reach back teeth.

Charging after two weeks



The AirFloss lasts two weeks between charging.

Microburst technology



Microburst technology applies a quick burst of pressurized air and micro-water droplets to clean deep between teeth where a toothbrush can't reach

Up to 99.9% Plaque Removal



AirFloss uses microdroplets of air and water or mouthwash to remove up to 99% more plaque in-between teeth*

Simple one-button use



Simple one-button use for easy cleaning in 60 seconds.

Slim, angled nozzle



Slim, angled nozzle for easy to even the hard-to-reach areas

Slim, easy to hold handle



Easy to manoeuvre around the mouth.

Specifications

Ease of use

Handle: Slim ergonomic design Battery indicator: Light in power button blinks to indicate when to recharge

Technical specifications

Battery: Rechargeable Charging time: 24 hour(s) Operating time (full to empty): 2 weeks or 14 full-mouth uses Operating time: Replace nozzle every 6 months

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty

© 2019 Koninklijke Philips N.V. All Rights reserved.

their respective owners.

Specifications are subject to change

Issue date 2019-07-23 Version: 4.0.1



