



AirFloss Pro/Ultra - Interdental cleaner

- 1 nozzle

HX8431/02

Healthier gums in 2 weeks*

Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro/Ultra**** is the easiest way to effectively clean between teeth. AirFloss Pro/Ultra can be used with mouthwash or water and is clinically proven as effective as floss for gum health.**

Benefits

Improves oral health

- Up to 99.9% plaque removal***
- Clinically proven as effective as floss for gum health**
- Helps prevent cavities between teeth

Superior technology

- Air and micro-droplet technology
- Triple burst customizable settings

- High performance nozzle

The easiest way to effectively clean between teeth

- Total daily confidence of a fresh and clean mouth
- An easy way to start a healthy habit
- Takes just 60 seconds: Point, Press, Clean!

PHILIPS


Features

Up to 99.9% plaque removal



AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.***

Micro-droplet technology



Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.** Helps improve gum health in as little as two weeks.

Fresh & clean with mouthwash



Fill the reservoir on the handle with either mouthwash or water, then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your liking, with single, double or triple bursts for each press of the activation button.

High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss Pro/Ultra is an easy way to clean deeper between teeth, helping to form a healthy habit.

Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple), and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

Specifications

[field not found 'disclaimer_text']

Items included

AirFloss Pro/Ultra handle	1
AirFloss Pro/Ultra nozzle	1
Charger	1

Cleaning performance

Clinically proven results	As effective as floss for gum health improvement for reducing gum bleeding for plaque reduction
Performance	Removes up to 99.9% plaque**
Health benefits	Cleans entire mouth in 60 sec Improve gum health in 2 weeks*
Customizable settings	Auto-burst or manual burst Single/double/triple burst***
For best results	Change nozzle every 6 months

Ease of use

Handle	Slim ergonomic design
Battery indicator	Handle LEDs show charge level
Battery Life	1-3 weeks****
Nozzle system	Easily snaps on and off

Technical specifications

Battery	Rechargeable
Battery type	Lithium ION
Charge time	24 hours to full charge

Power

Voltage	Multi-voltage charger
---------	-----------------------

Service

Warranty	2-year limited warranty
----------	-------------------------

Design and finishing

Color	Pink
-------	------

Packaging dimensions

Height	23.50 cm
Width	13.60 cm
Depth	8.40 cm
Nett weight	0.33 kg
Gross weight	0.46 kg
EAN	08710103866930
Number of products included	1
Country of origin	MX
Harmonised Systems Code	850980

Outer Carton

Length	26.80 cm
Width	24.91 cm
Height	15.70 cm
Gross weight	1.58 kg
EAN	18710103866937
Number of consumer packagings	3

** when used in conjunction with a manual toothbrush and anti-microbial mouth wash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under Support tab for further details.

*** From the treated areas; In a lab study, actual in-mouth results may vary

**** AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.

***** depending on burst setting used



data subject to change
2019, September 2

Version: 15.1.1
EAN: 08710103866930

© 2019 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com