

PHILIPS

sonicare

Interdental -
Rechargeable

AirFloss

w/ 3 Nozzles

Travel case

1 Nozzle holder



HX8283/08



Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard to reach areas.

Achieve a 30-second clean

- Cleans your entire mouth in just 30 seconds
- Point. Press. Clean!

Helps build a healthy habit

- An easy way to improve cleaning between teeth

Improves oral health

- Removes more plaque between teeth
- Helps improve gum health in just two weeks
- Helps prevent cavities between teeth

Easy to use

- Air and micro-droplet technology
- Comfortable ergonomic handle
- Guidance tip ensures correct placement



asimpleswitch.com


Highlights

An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

Removes more plaque between teeth



Removes more plaque between teeth

Air and micro-droplet tech



Air and micro-droplet technology

Guidance tip



The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

30 Second Clean



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from space to space. It will automatically burst every second.

Easy-fill reservoir



Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

Better gum health



Sonicare AirFloss is clinically proven to improve gum health in just two weeks.

Helps prevent cavities



By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

Ergonomic handle

Designed to be easy to maneuver in tight spaces, AirFloss has a comfortable ergonomic handle and a non-slip grip.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Ease of use

Handle: Slim ergonomic design

Battery indicator: Flashing indicates low battery

Battery Life: Up to 2 weeks

Nozzle attachment: Easily snaps on and off

Technical specifications

Battery: Rechargeable

Operating time: Replace nozzle every 6 months

Battery type: Lithium ION

Charge time: 24 hours to full charge

Power

Voltage: Multi-voltage charger

Service

Warranty: 2-year limited warranty

Items included

AirFloss Handle: 1

AirFloss Nozzle: 1, 3

Charger: 1 multi-voltage charger

Nozzle holder: 1

Cleaning performance

Performance: Removes up to 5x more plaque*

Health benefits: Improves gum health in only two weeks

Auto-burst feature: Cleans entire mouth in 30 sec

Nozzles: Best results change every 6mos

Design specifications

Color(s): White with green accents

