

# PHILIPS

## sonicare

Rechargeable powered  
dental flosser – pink nozzle

2 pink nozzles

237 ml mouthwash



HX8432/12



## Healthier gums in 2 weeks, guaranteed\*

### Designed for inconsistent flossers

For those who don't floss consistently, AirFloss Pro flosser is the easiest way to effectively clean between teeth. AirFloss Pro can be used with mouthwash or water and is clinically proven to be as effective as floss for gum health.\*\*

#### Improves oral health

- Up to 99.9% plaque removal\*\*\*
- Clinically proven to be as effective as floss for gum health\*\*
- Helps prevent cavities between teeth

#### Superior technology

- Air and micro-droplet technology
- High-performance nozzle
- Triple-burst customisable settings

#### The easiest way to effectively clean between teeth

- Total daily confidence with a fresh and clean mouth
- Takes just 60 seconds: Point, Press, Clean!
- An easy way to start a healthy habit

# Highlights

## Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.\*\*\*

## Micro-droplet technology




Our clinically proven results are possible thanks to our unique technology, which combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

## Easily start a healthy habit



Interproximal cleaning is very important to overall oral health. AirFloss Pro/Ultra is an easy way to clean deeper between teeth, helping to form a healthy habit.

## Better gum health



Philips Sonicare AirFloss Pro/Ultra nozzle is clinically proven to improve gum health as much as floss.\*\* Helps improve gum health in as little as two weeks.

## Helps prevent cavities




By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between your teeth.

## High-performance nozzle



The new AirFloss Pro/Ultra high-performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

## Triple-burst customisation



AirFloss Pro/Ultra's powerful spray can be customised to your liking, with single, double or triple bursts for each press of the activation button.

## Cleans in just 60 seconds



With the AirFloss Pro/Ultra, cleaning your entire mouth takes less than 60 seconds a day. Simply select your burst frequency (single, double or triple) and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

## Fresh and clean with mouthwash



Fill the reservoir on the handle with either mouthwash or water, and then point and press. Use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

# Specifications

## Items included

AirFloss Pro/Ultra handle: 1  
AirFloss Pro/Ultra nozzle: 2  
Mouthwash: 237 ml BreathRx  
Charger: 1

## Design and finishing

Colour: Pink

## Ingredients

**Mouth Rinse Ingredients:** Aqua, Sorbitol, Propylene Glycol, PEG-40, Hydrogenated Castor Oil, Xylitol, Poloxamer 407, Aroma, Zinc Gluconate, Cocamidopropyl Betaine, Cetylpyridium Chloride, Sodium Saccharin, Citric Acid, Limonene, BHT, CI42090

## Cleaning performance

**Clinically proven results:** As effective as floss, for gum-health improvement, for reducing gum bleeding, for plaque reduction

**Performance:** Removes up to 99.9% plaque\*\*

**Health benefits:** Cleans entire mouth in 60 secs, Improve gum health in 2 weeks\*

**Customisable settings:** Auto-burst or manual burst, Single/double/triple burst\*\*\*

**For best results:** Change nozzle every 6 months

## Ease of use

**Handle:** Slim ergonomic design

**Battery indicator:** Handle LEDs show charge level

**Battery Life:** Up to 2 weeks

**Nozzle attachment:** Easily snaps on and off

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Charge time:** 24 hours to full charge

## Power

**Voltage:** Multi-voltage charger

## Service

**Warranty:** 2-year limited warranty



\* or your money back

\*\* when used in conjunction with a manual toothbrush and anti-microbial mouthwash in patients with mild to moderate gingivitis; AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. Please see Q&A under the Support tab for further details.

\*\*\* From the treated areas; In a lab study, actual in-mouth results may vary

\*\*\*\* depending on burst setting used

\*\*\*\*\* AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.