

# PHILIPS

## sonicare

AirFloss Pro/Ultra -  
Dispense

AirFloss Ultra

2 nozzles

2 oz. BreathRx mouth rinse



HX8382/13



## A triple burst of better gingival health

### Efficacy meets compliance for infrequent flossers

Philips Sonicare AirFloss Pro/Ultra\*\*\*\* is clinically proven as effective as floss in improving gingival health.\* In fact, up to 97% of inconsistent flossers who tried AirFloss Pro/Ultra showed improved gum health.\*

#### Improves patients' oral health

- Up to 99.9% plaque removal\*\*\*
- Clinically proven as effective as floss for gum health\*\*
- Helps prevent cavities between teeth

#### Superior technology

- Air and micro-droplet technology
- High performance nozzle
- Triple burst customizable settings

#### The easiest way to effectively clean between teeth


- An easy way to improve cleaning between teeth
- Takes just 60 seconds: Point, Press, Clean!
- Fresh and clean with mouthwash
- An easy way to start a healthy habit

#### Superior results for patients' oral care needs

- Exclusive bonus items only available to dental professionals


# Highlights

## Up to 99.9% plaque removal




AirFloss Pro/Ultra removes up to 99.9% of plaque from treated areas.\*\*\*

## Micro-droplet technology




Our clinically-proven results are possible from our unique technology that combines air and mouthwash or water to powerfully yet gently clean between teeth and along the gumline.

## An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

## Better gum health



Philips Sonicare AirFloss Pro/Ultra is clinically proven to improve gum health as much as floss.\*\* Helps improve gum health in as little as two weeks.

## Cleans in just 60 seconds



With the AirFloss Pro/Ultra, patients can clean their entire mouth in less than 60 seconds a day. Patients simply select their desired burst frequency (single, double or triple) and hold down the activation button for continuous automatic bursts or press and release for manual burst mode.

## Fresh and clean with mouthwash




Patients can use with mouthwash for the ultimate fresh experience and anti-microbial benefits.

## Helps prevent cavities between



By gently bursting away plaque that brushing missed, Philips Sonicare AirFloss Pro/Ultra helps prevent cavities from forming in the spaces between patients' teeth.

## High performance nozzle



The new AirFloss Pro/Ultra high performance nozzle amplifies the power of our air and micro-droplet technology to be more effective and efficient than ever.

## Triple burst customization



AirFloss Pro/Ultra's powerful spray can be customized to your patient's liking with single, double, or triple bursts for each press of the activation button.

# Specifications

## Items included

AirFloss Pro/Ultra handle: 1  
AirFloss Pro/Ultra nozzle: 2  
Antibacterial mouth rinse: 2oz  
Nozzle holder: 1  
Travel bag: 1  
Charger: 1

## Design and finishing

Color: Gray

## Cleaning performance

**Clinically proven results:** As effective as floss, for gum health improvement, for reducing gum bleeding, for plaque reduction

**Performance:** Removes up to 99.9% plaque\*\*

**Health benefits:** Improve gum health in 2 weeks\*, Cleans entire mouth in 60 sec

**Customizable settings:** Auto-burst or manual burst, Single/double/triple burst\*\*\*

**For best results:** Change nozzle every 6 months

## Ease of use

**Handle:** Slim ergonomic design

**Battery indicator:** Handle LEDs show charge level

**Battery Life:** 1-3 weeks\*\*\*\*\*

**Nozzle attachment:** Easily snaps on and off

## Technical specifications

**Battery:** Rechargeable

**Battery type:** Lithium ION

**Charge time:** 24 hours to full charge

## Ingredients

**Antibacterial mouth rinse:** 0.075% CPC, Zytex

## Power

**Voltage:** Multi-voltage charger

## Service

**Warranty:** 2-year limited warranty



\* When used in conjunction with a manual toothbrush and antimicrobial rinse in patients with mild to moderate gingivitis. AirFloss is designed to help inconsistent flossers develop a healthy daily interdental cleaning routine. For more information, please visit [www.philips.com/airfloss/faq](http://www.philips.com/airfloss/faq).

\*\*\* \*AirFloss Ultra and Pro are the same product but may be named differently depending on the country and channel.  
\*\*\*\*\* depending on burst setting used