

# PHILIPS

## sonicare

Interdental -  
Rechargeable

AirFloss

w/ 3 Nozzles



HX8243



## Don't floss? Then AirFloss.

Removes plaque where brushing can't

For those who don't floss consistently, AirFloss is an easy way to start cleaning between teeth. AirFloss can be used with mouthwash or water and features unique air and micro-droplet technology to help remove plaque in hard to reach areas.

### Achieve a 30-second clean

- Cleans your entire mouth in just 30 seconds
- Point. Press. Clean!

### Helps build a healthy habit

- An easy way to improve cleaning between teeth

### Improves oral health

- Removes more plaque between teeth
- Helps improve gum health in just two weeks
- Helps prevent cavities between teeth

### Easy to use


- Air and micro-droplet technology
- Comfortable ergonomic handle
- Guidance tip ensures correct placement



[asimpleswitch.com](http://asimpleswitch.com)


# Highlights

## An easy way to start



Interproximal cleaning is very important to overall oral health. AirFloss is an easy way to clean deeper between teeth, helping to form a healthy habit. After three months of using AirFloss, 96% of inconsistent flossers continued to use AirFloss four or more days per week.

## Removes more plaque between teeth



Removes more plaque between teeth

## Air and micro-droplet tech



Air and micro-droplet technology

## Guidance tip



The slim angled nozzle and guidance tip make finding the correct placement easy. Just slide along the tooth's surface near the gum line until the tip fits into the groove between the teeth.

## 30 Second Clean



With the AirFloss, cleaning your entire mouth takes just 30 seconds. Simply press the button once for a single burst, then move and repeat. To use the new auto-burst feature, hold down the button and move the guidance tip from space to space. It will automatically burst every second.

## Easy-fill reservoir



Fill the reservoir on the handle with either mouthwash or water for a refreshing clean, then point and shoot. The reservoir holds enough liquid for up to two uses.

## Better gum health



Sonicare AirFloss is clinically proven to improve gum health in just two weeks.

## Helps prevent cavities



By gently bursting away plaque that brushing missed, Sonicare AirFloss helps prevent cavities from forming in the spaces between your teeth.

## Ergonomic handle

Designed to be easy to maneuver in tight spaces, AirFloss has a comfortable ergonomic handle and a non-slip grip.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Ease of use

**Handle:** Slim ergonomic design

**Battery indicator:** Flashing indicates low battery

**Battery Life:** Up to 2 weeks

**Nozzle attachment:** Easily snaps on and off

## Technical specifications

**Battery:** Rechargeable

**Operating time:** Replace nozzle every 6 months

**Battery type:** Lithium ION

**Charge time:** 24 hours to full charge

## Power

**Voltage:** Multi-voltage charger

## Service

**Warranty:** 2-year limited warranty

## Items included

**AirFloss Handle:** 1

**AirFloss Nozzle:** 1, 3

**Charger:** 1 multi-voltage charger

**Nozzle holder:** 1

## Cleaning performance

**Performance:** Removes up to 5x more plaque\*

**Health benefits:** Improves gum health in only two weeks

**Auto-burst feature:** Cleans entire mouth in 30 sec

**Nozzles:** Best results change every 6mos

## Design specifications

**Color(s):** White with green accents

