



Philips Sonicare Power
Flosser 7000
Oral Irrigator

Quad Stream technology

Whisper quiet design
Pulse waves guide as you go



HX3904/32

Removes up to 99.9% of plaque in treated areas*

Flossing reinvented for superior gingival health

With distinctive features that give patients a new level of clean and encourage lasting compliance, the new Philips Sonicare Power Flosser is the next wave in flossing.

New level clean

- Quad Stream for a thorough clean
- Pulse Wave technology guides from tooth to tooth
- Whisper quiet for great results without disturbing others
- 4 flossing modes, 10 intensities
- Easy magnetic dock with auto shut-off
- F1 Standard nozzle for a regular clean
- F3 Quad Stream nozzle for maximum coverage
- Easy to fill 600ml reservoir
- Compact nozzle case


PHILIPS
sonicare

Highlights

For a thorough clean

Unique cross-shaped nozzle creates four wide streams for a thorough clean interproximally and along the gingival margin.

Guides from tooth to tooth



Brief pauses in the cleaning pulsations prompt your patients to glide the nozzle to the next interproximal space to ensure a whole-mouth clean.

Whisper quiet design



Engineered for a discreet clean, so your patients get great results without disturbing others.


4 modes, 10 intensities



Your patients will enjoy a clean that's tailored to their needs. Clean mode uses a continuous


flow of water for a great everyday clean; Deep Clean uses distinct pulses for a more thorough clean; Sensitive mode is extra gentle for sensitive gums; Gum care mode soothes and stimulates the gumline. Your patients can adjust the intensity for their comfort using the LED light display as their guide.

F3 Quad Stream nozzle



For maximum coverage and an extra thorough clean, use the Quad Stream nozzle to concentrate water into four wide streams between teeth and along the gum line. The soft rubber guide is gentle on gums and every nozzle simply clicks on and off.

F1 Standard nozzle



Standard nozzle for power cleaning between teeth to remove plaque and debris

Easy magnetic dock



When the patient is finished, they can simply return the handle to the easy magnetic dock. The handle automatically switches from pause to off.

Compact nozzle case

The patients can store extra nozzles hygienically in a compact and stylish nozzle case.

600ml reservoir



The 600ml reservoir with easy flip lid holds enough water for one cleaning. To refill simply lift out, fill, and click back in.

Specifications

Modes

- Clean: For exceptional everyday clean
- Deep Clean+: For an invigorating deep clean
- Sensitive: Gentle teeth and gum cleaning
- Gum Care: Gently massages gums
- Intensities: 10

Cleaning performance

- Performance: Removes up to 99.9% plaque*

Items included

- Power Flosser: 1
- F1 Standard nozzle: 1

- F3 Quad Stream nozzle: 1
- Nozzle storage case: 1

Ease of use

- Handle: Slim ergonomic design
- Nozzle attachment: Easily clicks on and off

Power

- Voltage: Multi-voltage charger

Service

- Warranty: 2-year limited warranty



Issue date 2020-07-20

Version: 4.0.1

© 2020 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Up to 99.9% plaque removal in treated areas in a lab study. Actual in-mouth results may vary.