

PHILIPS

sonicare

Sonic electric
toothbrush with app

ExpertClean 7400

Connected brushing made easy

Built-in pressure sensor

Smart brush head recognition

3 modes, 3 intensities



HX9685/03



Improving patient routines for exceptional check-u

Expert technology. Personalized guidance.

The Philips Sonicare ExpertClean focuses on how patients brush, helping them improve their oral care routine in between visits.

Proven to improve oral health

- Increased surface contact for eliminating plaque
- Better gum care for overall health

Designed around your patient

- Guides patients to better oral care habits
- 3 modes, 3 intensities
- Brush heads automatically select optimal settings
- Guides patients to use correct amount of pressure
- For the most effective cleaning, every day
- Alerts patients of low battery

Improves patients' oral health

- Removes 10x more plaque* for a deep clean
- Up to 7x healthier gums* in just 2 weeks


Highlights

Removes up to 10x more plaque*




The C3 Premium Plaque Control brush head delivers our deepest clean. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact** and removes up to 10x more plaque from hard-to-reach areas*.

Up to 7x healthier gums*



The G3 Premium Gum Care brush head improves patient gum health. Smaller size and targeted gum line bristles deliver gentle, yet effective, cleaning along the gumline, where gum disease starts. It's clinically proven to deliver up to 100% less gum inflammation* and up to 7x healthier gums in just two weeks.*

Progress Report



Sensors track and measure patients' brushing behaviors while the Philips Sonicare ExpertClean brush handle provides instant feedback. Over time, the Sonicare app retains and analyzes these habits and generates a Progress Report to help promote better brushing technique.

3 modes, 3 intensities




With three brushing modes (Clean, Gum Health, DeepClean+) and three intensity settings (Low, Medium, High), your patients can personalize their brushing experience for optimal comfort and a complete clean. See Specifications for more detailed information on modes.

Smart Brush Head Recognition



A microchip in the brush head syncs with the Philips Sonicare ExpertClean handle. If a patient clicks on a C3 Premium Plaque Control/Defense/Defence, G3 Premium Gum Care or W3 Premium White brush head, the toothbrush automatically selects the optimal brushing mode and intensity level. All your patients need to do is start brushing. It helps patients with having the best settings for the best performance and their needs

Pressure sensor



If patients apply too much pressure, the intuitive pressure sensor will gently pulsate the handle and give patients the real-time feedback they need to adopt a softer touch. In

a study of U.S. patients, 7 out of 10 aggressive brushers tested reduced their brushing pressure when using the pressure sensor.****

C3 Premium Plaque Control



C3 Premium Plaque Control is our deepest cleaning brush head yet. With soft flexible sides, the bristles contour to the shape of each tooth's surface to provide 4x more surface contact and remove up to 10x more plaque from hard-to-reach areas.*

G3 Premium Gum Care brush head



G3 Premium Gum Care is our best brush head for improving gum health. Its smaller size and targeted gumline bristles deliver gentle yet effective cleaning along the gumline, where gum disease starts. It's clinically proven to deliver up to 7x healthier gums* in just 2 weeks.

Specifications

Modes

3 intensities: Low, Medium, High

Clean: For exceptional everyday clean

Gum Health: Special attention to molars

Deep Clean+: For an invigorating deep clean

Items included

Handle: 1 ExpertClean

Brush heads: 1 C3 Premium Plaque Control, 1 G3 Premium Gum Care, 1 W3 Premium White

Travel case: 1

Charger: 1

BrushSync Mode Pairing

C3 Premium Plaque Control: Pairs with Clean mode

G3 Premium Gum Care brush head: Pairs with Gum Health mode

Connectivity

Bluetooth® wireless technology: Connected brushing app

Compatibility

Android compatibility: Android phones, Bluetooth 4.0 enabled tablets

iOS compatibility: iPad 3rd Gen or higher, iPhone 4S or higher, with iOS7 or higher, with iOS7 operational system

Design and finishing

Color: White and Light Blue

Cleaning performance

Performance: Removes up to 10x more plaque*

Health benefits: Up to 7x healthier gums*

Whitening benefits: Helps naturally remove stains

Pressure feedback: Vibrates handle to alert user

Timer: BrushPacer and SmarTimer

Smart sensor technology

Pressure sensor: Alerts when brushing too hard

Location sensor: Tracks & improves coverage

Scrubbing sensor: Guides to reduce scrubbing

BrushSync Replacement Reminder: Always know when to, replace brush heads

BrushSync Technology: Connects smart handle and, smart brush head

Ease of use

Handle compatibility: Easy click-on brush heads

Replacement reminder: To always ensure best results, reminder icon lights up

Handle: Slim ergonomic design

Technical specifications

Battery: Rechargeable

Battery type: Lithium ION

Operating time (full to empty): 14 days***

Power

Voltage: 110-220 V

Service

Warranty: 2-year limited warranty

