

# PHILIPS

Patient monitoring

Medical Tablet

Distributed by Philips

## Easy, efficient access to the information you need

The Philips Medical Tablet is built for healthcare environments to provide access to relevant information where you need it.



### Philips Medical Tablet

#### Main features

- 25.65 cm (10.1 inch) color LCD, 1280 x 800 HD, capacitive multi-touch panel
- MIL-STD-810G & IP65 certified, 1.2 m (4 ft) drop test onto concrete
- Connections include WLAN, Bluetooth, NFC and LAN
- Battery life up to 7 hours
- Operating system: Windows 10
- 2 MP/8 MP Camera

## System

Processor	Intel® Celeron® Processor N3350 (2 M cache, up to 2.4 GHz)
System Memory	Onboard LPDDR4 4 GB

## Storage

Solid State Drive	128 GB solid state drive
-------------------	--------------------------

## Panel

LCD Size	25.65 cm (10.1 inch)
Display Type	WXGA
Resolution	1280 x 800 (HD)
Luminance	300 cd/m <sup>2</sup>
Touch Type	10 points projective capacitive

## I/O

USB	2 x USB 3.0 type A
HDMI	1 x HDMI output
Other	1 x audio jack, 1 x micro SD, 1 x DC jack

## System control

Control Buttons	1 x power button
LED Indicator	1 x power

## Mechanical & environmental

Operating Temperature (with table top dock or wall mount dock)	0 – 40°C (32 – 104°F)
Storage Temperature	-20 – 60°C (-4 – 140°F)
Operating Humidity	15 – 95% relative humidity, non-condensing
Weight	1 kg (2.2 lb)
Construction	Rubber + plastic
IP Rating	Dust-tight, and protected against splashing and water jets (IP65), and cleaning and disinfection agents
Vibration Test	MIL-STD-810G Method 514.6 Procedure 1
Shock Test	MIL-STD-810G Method 516.6 Procedure 1
Dimension	281 x 193.3 x 19.95 mm
Drop	1.2 m (4 ft) drop test on concrete floor, MIL-STD-810G 516.5 Procedure VI

## Certifications information and standards

Safety certification	CE-LVD, CCC, BSMI
EMC certification	CE-EMC, VCCI, FCC, BSMI, CCC
RF certification	CE-RED, FCC/IC, NCC, SRRC, MIC
Energy Efficacy certification	ErP
Medical electrical standards	IEC 60601-1, IEC 60601-1-2

## Software support

OS Information	Windows 10
----------------	------------

## Power requirement

Power Connector	DC jack
Power Requirement	19 V
Battery	28.44 W Li-ion battery
Battery Operating Time	7 hours (depending on usage)
Adapter	Input: 100 – 240 V AC / 50 – 60 Hz Output: 60 W adaptor (19 V@3.15 A adaptor)

## In-box accessory

Accessory	19 V 60 W power adaptor
-----------	-------------------------

## Data collection

Camera	2 MP front camera 8 MP rear camera with auto focus
NFC	ISO/IEC 14443 A/B, 15693/18092, NFCIP-1, NFCIP-2, Support MIFARE and Felica

## Wireless connections

WLAN	802.11 a/b/g/n/ac
LAN	(via docking accessories)
Bluetooth	Bluetooth 4.1

## Ordering information

	Order number
Philips Medical Tablet	867365
Medical Tablet	T01 (mandatory)
Table Top Dock	D01 (optional)
Wall Mount Dock	D02 (optional)

**To learn more about the Philips Medical Tablet and how it could be deployed at your facility, please contact your local Philips sales representative.**

