

PHILIPS

AVENT

Large Massage
cushion accessory

Fits all comfort breast pump



SCF167/01



Right size, maximum comfort

Breast pump with wider funnel

This soft massage cushion has a wider funnel, specially designed for moms with larger nipples. It can be used with Philips AVENT Comfort breast pump range. See below images to help you decide if you need the larger breast pump cushion.

Comfortable pumping

- Wider funnel for comfortable fit
- Fits on all Comfort breast pumps

Check your fit for optimal comfort

- Right fit: The nipple sits comfortably in the cushion funnel
- Choose large size: Nipple touches the walls of the funnel



asimpleswitch.com


Highlights

Large breast pump cushion



You can find out whether you need a large breast pump cushion by checking whether your nipple fits comfortably in the cushion funnel. If your nipple touches or rubs against the walls of the cushion funnel, then you should choose the large breast pump cushion to ensure a more comfortable expressing experience.

Correct fit



If the nipple sits comfortably in the cushion funnel, you are using the correct size massage cushion.

Too small fit



If your nipple touches or rubs against the wall of the cushion funnel, then you should choose the large breast pump cushion with wider funnel to ensure a more comfortable pumping experience.

Fits all Natural breast pumps

The large size soft massage cushion is compatible with all Philips AVENT Comfort breast pumps.

Specifications

Country of origin

England

Material

BPA free*

What is included

Large size cushion: 1 pcs

Cushion

Material: Silicone, BPA free*

Compatibility

Compatible with Natural range

Ease of use

Cushion: Easy to clean, Dishwasher proof

Development stages

Stage: Pregnancy, 0 - 6 months

Functions

Soft massage cushion: Gentle stimulation

Replaceable part

Fits product types :: SCF330

Fits product types:: SCF332

Fits product types :: SCF334



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.