



#### Baby Care set

- All baby care essentials
- Complete set
- Boys and girls

SCH400/00

## My first baby care set

All baby care essentials in one complete set

With our Baby care set SCH400 you can give the best care to your little one. The compact set includes an accurate, comfortable & quickly measuring digital thermometer, soft tip nasal aspirator, finger tooth brush, nail and hair care set.

#### Benefits

##### 10 baby care essentials

- Digital thermometer with professional accuracy\*
- Nasal Aspirator with soft and flexible tip
- Finger tooth brush
- Nail and hair care set

##### Compact pouch perfect for when travelling

- Compact and structured set
- A perfect gift

**PHILIPS**

## Features

### Digital thermometer



Quickly, comfortably and accurately measure your baby's temperature. The digital thermometer has professional accuracy\* and the flexible tip provides extra comfort for you and your baby.

### Compact and structured

The structured set offers room for your additional baby care products, whilst it is compact and convenient to take with you. Ideal for travel and nursery.

### Nasal aspirator



The nasal aspirator helps to open your babies airways to help your baby feel and sleep better.

## Specifications

[field not found 'disclaimer\_text']

### Technical specifications

|                   |                               |
|-------------------|-------------------------------|
| Battery           | Battery 1,55VDC (LR41)        |
| Battery life      | > 100 operating hours         |
| Measurement range | 35°C to 42.9°C range          |
| Accuracy          | ±0.1°C for 35°C to 42°C range |
| Compliance        | Conforms to MDD93/42/EEC      |

### Logistic data

|                              |                     |
|------------------------------|---------------------|
| F-box dimensions (W x H x D) | 117 x 174.5 x 51 mm |
|------------------------------|---------------------|

### Development stages

|       |            |
|-------|------------|
| Stage | All stages |
|-------|------------|

### Gender

|     |              |
|-----|--------------|
| For | Boy and girl |
|-----|--------------|

### Packaging dimensions

|        |          |
|--------|----------|
| Height | 17.70 cm |
|--------|----------|

### Finger tooth brush

Softly massage your baby's gums or teeth.

### Nail and hair care set



Complete nail care set with rounded-tip scissors, nail clipper and three emery boards. And a hair care set with rounded-tip comb and soft hair brush.

### Gift

A perfect gift for all (first time) parents.

|                             |                |
|-----------------------------|----------------|
| Width                       | 11.80 cm       |
| Depth                       | 5.20 cm        |
| Nett weight                 | 0.22 kg        |
| Gross weight                | 0.26 kg        |
| EAN                         | 08710103653783 |
| Number of products included | 1              |
| Country of origin           | CN             |
| Harmonised Systems Code     | 902519         |

### Outer Carton

|                               |                |
|-------------------------------|----------------|
| Length                        | 24.90 cm       |
| Width                         | 17.20 cm       |
| Height                        | 18.90 cm       |
| Gross weight                  | 1.72 kg        |
| EAN                           | 18710103832291 |
| Number of consumer packagings | 6              |

\*\* ±0.1°C between 35°C and 42 °C at room temperature 22°C



data subject to change  
2020, December 22

Version: 4.2.1  
EAN: 08710103653783

© 2020 Koninklijke Philips N.V.  
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

[www.philips.com](http://www.philips.com)