

PHILIPS

AVENT

Microwave steam sterilizer


Kills 99.9% of harmful germs

Sterilizes in 2 minutes

Fits 4 Philips Avent bottles

Fits most microwaves

SCF281/03



Fast microwave steam sterilizer

For use at home or away

The Philips Avent microwave steam sterilizer is a great choice for fast, effective sterilizing at home or away. Up to 4 Philips Avent bottles or products can be sterilized at once, with 99.9% of germs killed in just 2 minutes.

Fast and effective microwave steam sterilizing

- Microwave steam sterilizer ready in just 2 minutes
- Kills 99.9% of germs and bacteria
- Contents stay sterile for up to 24 hours

Fits various bottle types and other products

- Sterilizes standard neck and wide neck baby bottles
- Sterilizes pumps, soothers, cutlery as well as baby bottles

Compact, lightweight design

- Lightweight design for sterile baby bottles on the go
- Fits most microwaves on the market
- Stay cool safety clips to keep lid secure



asimpleswitch.com

Highlights

Sterilizes in just 2 minutes



The microwave steam sterilizer can be used to sterilize baby bottles and other products in just 2 minutes, killing 99.9% of germs and bacteria. The exact length of the cycle depends on the wattage of your microwave, of course. 2 minutes at 1200-1850W, 4 minutes at 850-1100W, 6 minutes at 500-800W.

Kills 99.9% of germs



The microwave steam sterilizer is proven to kill 99.9% of germs and bacteria.

Effective sterilization



Contents remain sterile for up to 24 hours if the lid is unopened.

Works across bottle types



The sterilizer sterilizes standard neck and wide neck baby bottles. Up to 4 Philips Avent baby bottles can be sterilized at once.

Use for pumps, soothers too

The sterilizer can be used to sterilize breast pumps, soothers, cutlery and other products as well as baby bottles.

Lightweight design for travel



The sterilizer's lightweight and compact design makes it easy to take along with you wherever you go - holidays, visits to relatives etc. Now you can always have sterilized baby bottles and other products within easy reach.

Fits most microwaves



The microwave steam sterilizer has been designed to fit practically all microwaves on the market, keeping baby bottles and other products sterile wherever you are.

Stay cool safety clips



The microwave steam sterilizer has stay cool safety clips to keep the lid secure. This ensures you can remove the sterilizer safely from the microwave, keeping baby bottles and other products sterile and preventing the possibility of burning.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas - Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Technical specifications

Sterilization time: 2 min at 1200-1850W, 4 min at 850-1100W, 6 min at 500-800W
Water capacity: 200ml

Weight and dimensions

Weight: 740 g
Dimensions: 166 (H), 280 (W), 280 (L) mm

Country of origin

Made in: Poland

What is included

Microwave steam sterilizer: 1 pcs
Tongs: 1 pcs

Development stages

Stages: 0 - 6 months

