

PHILIPS

InfraCare infrared
lamp

300W

Zone treatment area

Digital timer



HP3631



Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 300 Watt infrared lamp is very comfortable and the infrared warmth penetrates deep into the skin, stimulates blood circulation and gives relief for "zone" areas of 40x30 cm.

Effective pain relief

- Deep penetrating Infrared warmth

Zone treatment

- 300 Watt infrared halogen lamp

Easy to use

- Adjustable angle

Treatment control

- Digital timer to set personal treatment time

Highlights

300 Watt infrared halogen lamp



The Philips InfraCare HP 3631 has been optimised for effective relief of pains in areas covering 40 x 30 cm, such as both shoulder and neck, the lower back, a thigh, an arm. The InfraCare features innovative infrared halogen lamp technology. The special optics, filter and the powerful halogen lamp of 300 Watt have been optimised, enabling zone treatment (40 x 30cm), and a much more comfortable warmth, as the infrared warmth is more equally divided over the treatment area.

Adjustable angle



The infrared lamp can easily be positioned upto 40 degrees backwards.

Digital timer to set personal treatment time



Treatment time can easily be controlled by setting the advised or preferred number of minutes. The digital timer automatically counts down and gives auditative feedback at the end of the session.

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

Specifications

Effective pain relief

Deep penetrating warmth

Zone treatment

Treatment area: 40x30 cm

Infrared halogen lamp: 300 W

Easy to use

Adjustability: 0-40 degrees (backwards)

Soft-touch handgrip: for easy transport and positioning

On/off switch

Technical specifications

Power: 300 W

Voltage: 230 V

Frequency: 50/60 Hz

Lifetime of lamps: 500 hour(s)

Battery: LR-54-Alkaline, 1,5 V (lifetime 200 actions)

Cord length: 200 cm

Insulation: Class II (double isolation)

Weight and dimensions

Product dimensions: 20x30.9x17.2 (WxHxD) cm

Product weight: 1.4 kg

F-box dimensions: 28x38.6x24.5 (WxHxD) cm

F-box weight: 1.8 kg

No. of F-boxes in A-box: 3 pcs

A-box dimensions: 29x40.4x74.5 (WxHxD) cm

A-box weight: 5.4 kg

Qnt. on Euro pallet: 60 pcs

Logistic data

CTV code: 884363101000

Country of origin: Germany

Safety

IEC certified: Complies to IEC 60601 and 60335

Treatment control

Digital timer: with auditative feedback

Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

