

# PHILIPS

## AVENT

Anti-colic with  
AirFree™ vent

3 Bottles

9oz/260ml

Slow flow teat

1m+



SCF403/35



## Designed to reduce colic, wind and reflux\*

### Locks milk in, keeps air out

The AirFree™ vent is designed to help your baby swallow less air. The teat stays full of milk, even when horizontal, for upright feeding. Reducing the air a baby ingests can help ease common feeding issues like colic, gas and reflux.

#### Designed for less air intake

- Locks milk in, keeps air out, for easy upright feeding
- Teat stays full of milk, not air

#### Clinically proven to reduce colic and fussing\*

- Anti-colic system proven to reduce colic and fussing\*

#### Easy to use -- and clean!

- Easy to clean & assemble, as AirFree™ vent is a single piece
- Simple to use
- Teat shape designed for a secure latch-on

#### Additional benefits

- Use your bottle with or without the AirFree™ vent
- Different teat flow rates available
- Mix and match with other Philips Avent products



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Locks milk in, keeps air out



The teat stays full of milk, even when the bottle is horizontal, so your baby can drink in a more natural, upright position. This can help reduce reflux, aid digestion, and make feeding time more comfortable for you and your little one.

## Designed for less air intake



Our unique AirFree™ vent draws air away from the teat, so your baby swallows less air as they drink. This can help reduce common feeding issues such as colic, reflux and wind.

## Reduces colic and fussing\*



Clinical studies have shown that the Philips Avent bottle reduces colic and fussing\*. In addition, the valve in the teat prevents vacuum build-up. This means that the air stays in the bottle, away from your baby's tummy.

## Easy to clean and assemble



The AirFree™ vent can be easily assembled with any Philips Avent Anti-colic and Classic+ bottle. Also cleaning is a breeze, thanks to the bottle's wide-neck, and there only being a few parts to clean.

## Simple to use



To start, fill the bottle with milk, place the AirFree™ vent over the rim, and screw on the teat. Tilt the bottle downwards once to fill the teat, and you're ready to go. When feeding, make sure the AirFree™ vent stays on top. The teat will stay full of milk, even when it's horizontal, for easier upright feeding.

## Designed for a secure latch




The shape of the teat allows for a secure latch-on, and the ribbed texture helps prevent teat collapse, for uninterrupted, comfortable feeding.

## Freedom of choice



You can use your Philips Avent Anti-colic bottle with or without the AirFree™ vent. The AirFree™ vent is compatible with our 4oz/125ml, 9oz/260ml and 11oz/330ml Philips Avent Anti-colic and Classic+ bottles.

## Choosing the right teat



As your baby grows, their drinking needs change. Choose a teat with a flow rate that's right for their development. The Philips Avent Anti-colic and Classic+ bottles can be used with or without the AirFree™ vent—but they should only be used with Philips Avent Anti-colic and Classic+ teats.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Design

**Bottle design:** Ergonomic shape, Wide neck

## Development stages

**Stage:** 0-12 months

## Ease of use

**Bottle use:** Easy to assemble, Easy to clean, Easy to hold, Dishwasher & microwave safe

## Functions

**Ease of use:** Easy to clean and assemble, Few parts for easy assembly

**Teat:** Easy latch on, Ribbed texture prevents teat collapse, Clinically proven Anti-colic system\*

**Compatibility:** Full Anti-colic and Classic+ bottle compatibility

## Material

**Bottle:** Polypropylene, BPA free

**Teat:** Silicone, BPA free

## What is included

**Baby bottle:** 3 pcs

**AirFree™ vent:** 3 pcs



\* The AirFree™ vent is designed to help your baby swallow less air. The teat stays full of milk, even when horizontal, for upright feeding. Reducing the air a baby ingests can help ease common feeding issues like colic, wind and reflux.

\* 80% of moms agreed that "my baby experienced fewer feeding issues" in a home placement test with 144 mothers in the US in 2017

\* At 2 weeks of age, babies fed with a Philips Avent bottle showed a trend to less colic compared to a conventional bottle and a significant reduction in fussing at night compared to babies fed with another leading bottle.

\* What colic is, and how it affects babies? Colic is caused in part by swallowing air while feeding, which creates discomfort in a baby's digestive system. Symptoms include crying, fussing, wind and spit-up.