

# PHILIPS

# AVENT

## Easy Comfort Single Electric Breast Pump

Includes one 120 ml/4 oz bottle



SCF301/03



## Convenience and Comfort for effective expression

Choose between 4 stimulation modes and 4 pumping settings to make every expression session comfortable and personalised. The lightweight, compact design also means expressing on-the-go is a breeze.

### Multiple hassle-free settings for personalised use

- 4 stimulation modes x 4 pumping modes for optimal expressing
- Hassle-free, set and forget settings

### Designed to gently stimulate milk flow

- Soft massage cushion designed to help stimulate milk flow

### Lightweight and compact design

- Easy expressing, at home or on the go
- Battery and electricity operated

### Rhythmic Wave Technology

- Designed to mimic babies' suckling rhythm

### Discreet and hygienic expressing

- Discreet sound and backlit display for discreet use, day/night
- Hygienic expressing



[asimpleswitch.com](http://asimpleswitch.com)

# Highlights

## Personalised pumping



It's important to be comfortable when expressing and comfort means something different for every mum. With this breast pump, you can choose between 4 stimulation modes and 4 pumping settings. This way, you can make each expressing session efficient and comfortable.

## Hassle-free settings



Mums have enough on their mind without having to recall their ideal setting - so we'll do that for you. Your preferred settings are stored, so at the press of a button you're ready to go. Need a break from pumping? Simply press the pause button, and to resume press play. It's that easy.


## Soft massage cushion

Our breast pumps were designed with your comfort in mind. The soft and velvety texture of the massage cushions feel warm and comfortable against your skin, which can help stimulate milk flow. And the flexible material is designed to mimic your baby's suckling action, which can help milk letdown.

## Lightweight and compact design

The lightweight and compact design means you can easily express at home or take it with you for expressing on the go.

## Dual Power



The product runs on electricity. For extra mobility you can also insert batteries to allow cordless use. You can insert four non-chargeable 1.5 V AA batteries or four rechargeable 1.2 V AA batteries

## Designed to mimic suckling



When babies breastfeed, they follow a rhythm of suck, swallow, breathe, with a ratio of 1:1:1. This breast pump features special Rhythmic Wave Technology that's designed to mimic a baby's suckling and swallowing pattern, pace and ratio (1:1), to quickly and efficiently generate milk flow. This means your expressing sessions feel natural and comfortable.

## Discreet expressing



Our handy backlit display means you can express discreetly at night, without having to turn on the lights. And the soft sound means less disturbance.

## Hygienic expressing



We know that every drop counts, that's why we make sure all breastmilk flows straight from breast to container, without any backflow into the tubes. All parts that come into contact with your breastmilk are made using BPA-free materials.



[asimpleswitch.com](http://asimpleswitch.com)

## Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas — Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

# Specifications

## Country of origin

China

## Material

**Breast pump:** BPA-free\* (food contact parts only)

**Bottle:** BPA-free\*, Polypropylene

**Teat:** BPA-free\*, Silicone

## What is included

**Breast pump body:** 1 pcs

**Base unit incl. tubing:** 1 pcs

**Soft cushion:** 1 pcs

**Natural bottle 125 ml/4 oz:** 1 pcs

**Travel cover:** 1 pcs

**Sealing disc for milk storage:** 1 pcs

**Breast pad sample packs:** 1 pcs

**Lanyard:** 1 pcs

## Design

**Breast pump design:** Compact design

**Bottle design:** Ergonomic shape, Wide neck

## Ease of use

**Breast pump use:** Intuitive assembly, Easy cleaning, Fully compatible range

## Development stages

**Stage:** 0–6 months

## Functions

**Rhythmic wave technology:** Mimics baby suckling pattern

