

PHILIPS

AVENT

Natural Twin electrical
breast pump

Includes a 4oz/120ml Bottle

4 expression settings



SCF303/01



Comfort & convenience for effective expression

Rhythmic wave technology like baby's suckling

Our pump is designed to give you closest natural comfort feeling with convenience for fast and effective milk expression

Effective expression through more comfort

- Rhythmic wave technology just like baby's suckling
- Alternating suction delivers even suction pace
- Soft massage gently stimulates milk flow
- Optimized suction level offers just right suction

Effective expression through more convenience

- Memory/ pause helps store preferred settings, take a break
- Hygienic, safe materials and BPA free.
- Long life motor gives you high quality and durability
- Backlit is easy to express at night, no intrusive sound
- Clear digital display. Select modes, comfort level and more
- Lightweight including handy bag for easy transport
- Natural latch on and easy to combine with breastfeeding



asimpleswitch.com


Highlights

Rhythmic wave technology



Our rhythmic wave technology creates an expression cycle just like the way a baby naturally draws and swallows milk at the breast for fast and effectively generating milk flow. Fast and shallow stimulating movements are followed by slow and deep ones with an even natural suckling pace to gently but effectively express more milk. Pumping time is reduced and the expressing experience feels more natural and comfortable.

Alternating suction mode



Our new twin electric breast pump uses our alternating suction rather than simultaneous suction which delivers even suckling pace in a natural way. This gives one breast time to rest, while the other expresses, helping to care for breasts during expression. Alternating suction is more gentle on the breasts, helping make expression more comfortable and relaxing.

Soft massage cushion

The pump's unique massage cushion has a soft velvety texture that gives a warm feel to the skin to comfortably and gently stimulate milk flow. During pumping, it carefully mimics your baby's suckling for stimulating fast milk let-down. The soft massage cushion is designed to be comfortable and caring for attachment to breasts.

Optimized suction level

Our new twin electric breast pump is designed to offer just the right suction level. Choose from four stimulation modes to stimulate faster let-knowwhile the four expression modes make it easy and comfortable to quickly express milk

Memory, time and pause/play




Our new twin electric breast pump is designed to be completely intuitive to use. If you need to take a break during expressing, simply press the pause button and the play button to resume. The automatic memory function also makes it simple to store preferred settings so you can always express in the way that works best for you. The memory function will always show the last intensity level that was used. You can also manually select if you prefer.

Safe and hygienic expression



All products that come in contact with breastmilk are made using BPA-free safe materials, completely free from any potentially harmful chemicals. The pump is also specially designed so that there is no contact between the milk and the tube.

Long life motor, dual power



Our breast pump provides high quality in material and design with long life motor which can support pumping with durable performance. Dual power with both 4AA batteries and electrical power support which can allow moms to express anywhere and anytime.

Backlit display, quiet motor



Backlit display is easy to express even at night with no more annoying or loud pumping sounds. Our new twin electric breast pump is designed to be quiet across all setting levels. Its low level noise production makes it possible for you to express discreetly at any time, anywhere, regardless at any setting you choose.



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.

Specifications

Country of origin

China

Material

Bottle: BPA free*, Polypropylene

Breast pump: BPA free*, Polypropylene

Nipple: BPA free*, Silicone

What is included

Breast pad sample packs: 2 Day pads

Natural bottle 4oz: 2 bottles

Sealing disc for milk storage: 2 pcs

Standard size cushion: 2 pcs

Base unit incl. tubing: 1 pcs

Breast pump body: 2 pcs

Spare diaphragm: 1 pcs

Travel cover: 2 pcs

Breast pump

Material: BPA free*

Design

Bottle design: Ergonomic shape, Wide neck

Breast pump design: Compact design

Ease of use

Breast pump use: Fully compatible range

Development stages

Stage: Pregnancy, 0 - 6 months

Functions

Rhythmic wave technology: close to a baby's natural suckling for more milk in less time

Alternating suction mode: expression more comfortable and relaxing

Soft massage cushion: gentle stimulation

Setting: 4 expression settings and 4 stimulation settings

