

PHILIPS

AVENT

Single electric breast pump

Soft massage cushion (19.5mm)

Includes bottle and nipple

Sealing disc for storage



SCF332/33



More comfort, more milk*. Anytime, anywhere

Single electric breast pump with massage cushion

Sit comfortably with no need to lean forward and let our soft massage cushion gently stimulate your milk flow. Our quiet pump can be used anytime and anywhere, also with batteries. It is easy to set up, personalise, use and clean.

Clinically proven less leaning forward*

- Designed to sit comfortably without leaning forward

Comes with a soft massage cushion

- The soft massage cushion gently stimulates your milk flow

Quiet pump for discreet expression, anywhere

- Quiet pump for discreet expression, anywhere

Effective expression with 3 customisable settings

- Features a gentle stimulation mode and 3 pumping settings

Easy to set up, express and clean

- Easy to set up, express and clean. Closed system.
- Includes our Natural bottle and teat for natural latch on

Highlights

Relaxing position



The breast pump has a unique design, so your milk flows directly from your breast into the bottle while you sit up straight. As you pump in more comfort without leaning forward your milk flows more easily. When you pump in comfort it reduces muscle tension and fatigue.

Soft massage cushion



Our soft and velvety massage cushion warms your skin for comfortable, gentle stimulation of milk flow. The iconic petal cushion mimics baby suckling to gently help stimulate let down. Made of medical silicone, it ensures soft contact with the breast. With each cycle, the pump mimics massage techniques like finger tapping and kneading. Massage: reduces swelling; improves blood circulation and stimulation and facilitates milk lymph and milk duct drainage, improve stimulation to breast tissue and increases milk volume.*

Quiet pump



Designed for discreet expression, this quiet pump will let you express anywhere. As the pump is compact, quiet and discreet to use, you can take it with you anywhere, allowing you to express milk at your own convenience and maintain your milk supply.

Simple settings to choose from



The pump automatically starts in gentle stimulation mode to get your milk flowing. Once your milk starts flowing, you can adapt the suction rhythm by choosing from 3 pumping settings to make milk flow most comfortable for you. Depending on your own personal comfort, you may wish to use a higher suction mode that provides deeper suction to help you express more breast milk.

Easy set up, express and clean



With only a few parts and intuitive design, the pump is easy to assemble and use so you can easily express breastmilk for your baby, anytime. The closed system is designed for clean expressing, preventing milk spillage in the tubes. All parts that are in contact with breast milk are dishwasher-safe and the durable pump is designed to last a long time.

Includes our Natural bottle



Our Natural bottle and nipple are included with the breast pump. The bottle's soft nipple closely resembles the breast. And the wide breast-shaped teat has a flexible spiral design and comfort petals to help with natural latch on. This makes it easy for you and your baby to combine breast and bottle feeding.

Specifications

Design

Bottle design: Ergonomic shape, Wide neck
Breast pump design: Compact design

Development stages

Stage: 0 - 6 months

Ease of use

Breast pump use: Easy cleaning, Fully compatible range, Intuitive assembly

Functions

No leaning forward: Sit in a comfortable position

Settings: 1 Stimulation mode, 3 Expression settings

Soft massage cushion: Gentle stimulation

Material

Bottle: BPA free*, Polypropylene

Nipple: BPA free*, Silicone

Breast pump: BPA free* (food contact parts only)

What is included

Breast pump body: 1 pcs

Base unit incl. tubing: 1 pcs

Standard size cushion(19.5 mm): 1 pcs

Natural bottle 4oz: 1 pcs

Travel cover: 1 pcs

Sealing disc for milk storage: 1 pcs

Disposable breast pads: 4 pcs

Country of origin

England

