

PHILIPS

AVENT

Manual breast
pump

with breast milk containers



SCF290/20



Express more. Quickly.

Enhanced efficiency from a manual breast pump

A manual breast pump offers you total control of the suction level and pumping rhythm. The unique Philips Avent Manual Breast Pump assures you enhanced efficiency, being clinically proven to express more milk than a hospital electric pump.

Stimulate and Express

- Soft massage cushions trigger natural let-down, like baby

Clinically proven effective

- Patented Let-down Massage Cushion

Flexible and portable

- Easy to assemble and use

Highlights

Soft massage cushions



Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

Patented Let-down Massage Cushion



Five, soft petal massagers of the Philips Avent breast pump gently flex in and out as you pump and work with the vacuum to imitate your baby's suckling.

Easy to assemble and use



Quiet and portable.

Specifications

What is included

Manual breast pump: 1 pcs
Breast milk container (125 ml/ 4 oz): 2 pcs
Newborn nipple travel pack: 1 pcs
Bottle stand/funnel cover: 2 pcs
Spare parts: 1 pcs
Extra soft newborn flow nipple: 1 pcs
Sealing disc for milk storage: 2 pcs

Country of origin

England

Development stages

Stage: 0 - 6 months

