

PHILIPS

 Product Picture

AVENT

Manual breast
pump

with VIA Storage System

SCF300/12

Express more. Store easily.

Express, store and feed in easy steps

The unique Philips Avent Manual Breast Pump SCF300/12 offers you total control of the suction level and pumping rhythm. It assures you enhanced efficiency, being clinically proven to express more milk than a hospital electric pump.

Benefits and features:

- Gently stimulates fast milk flow

Clinically proven effective

- Patented Let-down Massage Cushion

Fits all Philips Avent breast pumps and teats

- Use adaptor to express and feed

For safe storage and transport

- Leak-proof twist-on system

Highlights

Gently stimulates fast milk flow

As you pump the patented Philips Avent Let-down Massage Cushion flexes in and out, gently massaging the area around the nipple to help stimulate natural milk flow.

Patented Let-down Massage Cushion



Five, soft petal massagers of the Philips Avent breast pump gently flex in and out as you pump and work with the vacuum to imitate your baby's suckling.

Leak-proof twist-on system

For safe storage and transport

Adaptor to express and feed

Use the adaptor to express milk and feed with the same cup.

Specifications

Country of origin

England

Development stages

Stage: 0 - 6 months

Material

BPA free*

What is included

Milk storage cup lid: 5 pcs

Milk storage cup adaptor: 2 pcs

Extra soft Newborn Flow nipple: 2 pcs

Newborn nipple travel pack: 2 pcs

Spare parts: 1 pcs

Pre-sterilized VIA Cup (180 ml/ 6 oz): 5 pcs

Manual breast pump: 1 pcs

