

PHILIPS

 Product Picture

AVENT

Manual breast
pump

with disposable breast pads

SCF290/15

Express more. Quickly.

Enhanced efficiency from a manual breast pump

A manual breast pump offers you total control of the suction level and pumping rhythm. The unique Philips Avent Manual Breast Pump assures you enhanced efficiency, being clinically proven to express more milk than a hospital electric pump.

Benefits and features:

- Clinically proven
- Includes 20 Avent Ultra Comfort Disposable Breast Pads

Helps you breastfeed longer

- More milk, naturally

Clinically proven effective

- Patented Let-down Massage Cushion

Stimulate and express

- Soft massage cushions trigger natural let-down, like baby

Flexible and portable

- Easy to assemble and use

Highlights

Easy to assemble and use



Quiet and portable.

Clinically proven

Clinically proven to express more milk than a hospital grade double electric pump*.

Patented Let-down Massage Cushion



The five soft petal massagers of the Philips Avent breast pump gently flex in and out as you pump and work with the vacuum to imitate your baby's suckling.

Soft massage cushions



The patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and are intended to help stimulate a fast natural let down.

Includes 20 Ultra Comfort Disposable Breast Pads

Unique four layered construction for ultra absorbency, comfort and dryness.

More milk, naturally



Expression allows breast milk to be available, even when you can't be there

Specifications

Country of origin

England

What is included

- Bottle stand/funnel cover: 2 pcs
- Breast milk container (125 ml/4 oz): 2 pcs
- Disposable breast pad: 20 pcs
- Extra soft Newborn Flow teat: 2 pcs
- Newborn nipple travel pack: 1 pcs
- Spare parts: 1 pcs
- Manual breast pump: 1 pcs

Development stages

Stage: 0–6 months

