

PHILIPS

AVENT

Twin Expression Kit

SCF308/02



Inspired by nature

Express more. Store easily

Personal pump parts for use with the Philips Avent Twin electric breast pump.

Benefits and features:

- Ultra comfortable and effective
- Personal pump parts for use with the Twin electric breast pump

Helps you breastfeed longer

- More milk, naturally

Stimulate and express

- Soft massage cushions trigger natural let-down, like baby
- Gentle vacuum mimics baby's suckling for steady milk flow

Highlights

Gentle-draw vacuum

The gentle-draw vacuum mimics baby's suckling to create a steady milk flow that requires less pumping.

Patented Let-down Massage Cushion



Five soft petal massagers gently flex in and out, working in tandem with the reliable vacuum. This is intended to help stimulate a fast, natural milk flow.

More milk, naturally



Expression allows breast milk to be available, even when you can't be there

Personal pump parts/accessory

This Twin expression kit's funnels are only compatible with the Twin electric breast pump, due to different valves and a closed diaphragm. The Twin expression kit will not fit the manual or single electric breast pump

Soft massage cushions



The patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and are intended to help stimulate a fast natural let down.

Specifications

Country of origin

England

What is included

Milk storage cup lid: 5 pcs

Storage Cup (240 ml/8 oz): 5 pcs

Personal pump parts: 2 pcs

Development stages

Stage: 0–6 months

