

PHILIPS

InfraPhil infrared
lamp

150 W

 Product Picture

HP3616



Relieves muscular pains

Deep, focused infrared warmth

Deep, focused warmth from infrared lamp

Easy to use

- On/off switch
- Adjustable angle

Relieves muscle pain

- Deep penetrating Infrared warmth

Easy to store

- Cord storage

Specifications

Easy storage

Cord storage

Technical specifications

Voltage: 220/230 or 240/250 V V

Power: 150 W

Cord length: 180 m

Frequency: 50 Hz

Insulation: Class II (double isolation)

Type of lamps: PAR 38E, 150 W + prismatic rings for more focus

Lifetime of lamp: 750 sessions of 10 minutes

Consumer Trade Item

Height: 22.5

Width: 20.4

Length: 14

Net Weight: .756

Gross Weight: .885

GTIN: 08710103652427

Country of origin: HU

Harmonized system code: 901820

Outer Carton

Height: 24

Width: 28.8

Length: 62.4

Net Weight: 4.536

Gross Weight: 5.714

GTIN: 18710103652424

Weight and dimensions

Product dimensions: 22.0 x 13.0 x 18.5 (H x W x D) cm

Product weight: 1 kg

A-box dimensions: 24.0 x 28.8 x 62.4 (H x W x D) cm

F-box weight: 1 kg

A-box weight: 6.4 kg

Qnt. on Euro pallet: 192 pcs

No. of F-boxes in A-box: 6 pcs

F-box dimensions: 22.5 x 14.0 x 20.4 (H x W x D) cm

Logistic data

CTV code: 884361601000

Country of origin: Germany

Safety

IEC certified: Complies to IEC 60335

Easy to use

Adjustability: 0–40 degrees (backwards)

On/off switch

Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

Other

Height: 207

Width: 100

Length: 120

GTIN: 68710103652429

Highlights

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

Adjustable angle

The infrared lamp can easily be positioned up to 40 degrees backwards.

On/off switch

Easy on/off switch at the back of the infrared lamp, so the appliance does not have to be unplugged after each treatment.

Cord storage



The infrared lamp has a cord storage holder at the back of the appliance. Simply wrap the cord around it and it can easily be stored again.

