

Philips Avent Gift Set
Natural Beginnings

Classic



SCD236/00



Breastfeeding and sterilising essentials

A convenient feeding set

Includes Microwave Steriliser, Manual Breast Pump, Feeding Bottles, Disposable Breast Pads and Sealing Discs

Express II Microwave Steriliser

- Ultra fast, easy to use
- Convenient for travel
- Compact and lightweight

ISIS Manual Breast Pump

- Clinically proven
- Gently stimulates fast milk flow
- Patented Let-down Massage Cushion

Disposable Breast Pads

- Unique four layered construction for maximum dryness
- Soft with discreet, contoured shape

Airflex Feeding Bottle

- Works with baby's natural feeding action
- For healthy, active feeding

PHILIPS
AVENT

Highlights

Airflex Feeding Bottle

A clinical trial demonstrated that at 2 weeks of age, babies fed with the Avent Bottle experienced less colic than babies fed with a conventional bottle. (www.philips.com/Avent).

ISIS Manual Breast Pump

In a clinical study, ISIS was proven to be as effective as a hospital electric pump. ISIS has an award-winning design that imitates your baby's suckling. The result - more milk faster and more naturally than with any other manual breast pump.

Express II Microwave Steriliser

Sterilises 6 bottles in 2 minutes. Just add water, load and place in the microwave. 2 min at 1100-1850 Watt, 4 min at 850-1000 Watt, 6 min at 500-800 Watt.

Disposable breast pads



Unique four layered construction for maximum dryness: 1. Ultra-soft top layer with nipple indent — keeps breast dry at all times 2. Ultra-absorbent layer, cushioned for extra comfort 3. Ultra-absorbent core — draws in moisture to guard against embarrassment 4. Ultra breathable outer layer — helps avoid sore nipples. Anti-slip with adhesive tape — keeps pad in place

Convenient for travel

Convenient for travel

Compact and lightweight

Compact and lightweight

Soft with discreet, contoured shape

Philips Avent breast pads have a contoured shape that allows you to wear them discreetly under your clothing.

Gently stimulates fast milk flow



Gently stimulates fast milk flow

Patented Let-down Massage Cushion



Patented Let-down Massage Cushion

For healthy, active feeding

For healthy, active feeding

Specifications

Country of origin

- England

Development stages

- Stage: 0–6 months

What is included

- ISIS Manual Breast Pump: 1 pcs
- Disposable breast pad: 40 pcs
- Sealing disc for milk storage: 2 pcs

What is included

- Feeding Bottle (260 ml/9 oz): 2 pcs
- Feeding Bottle (125 ml/4 oz): 2 pcs
- Express II Microwave Steam Steriliser: 1 pcs



Issue date 2021-02-12

Version: 9.4.1

© 2021 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com