

# PHILIPS

InfraPhil infrared  
lamp

100 W

HP1511



## Relieves muscular pains

Deep, focused warmth

Philips infrared lamp gives pain relief for muscles and joints. Its 100 Watt infrared warmth penetrates deep into skin, stimulates blood circulation and gives relief for areas of 20x30 cm.

# Specifications

## Easy to use

**Adjustability:** 0–40 degrees (backwards)

## Weight and dimensions

**Product dimensions:** 20.0 x 11.7 x 16.7

(HxWxD) cm

**Product weight:** 0.5 kg

## Technical specifications

**Power:** 100 W

**Cord length:** 160 m

**Voltage:** 220/230 V

**Frequency:** 50 Hz

## Weight and dimensions

**A-box dimensions:** 27.0 x 24.0 x 55.0

(HxWxD) cm

**F-box weight:** 1 kg

**No. of F-boxes in A-box:** 6 pcs

**Qnt. on Euro pallet:** 252 pcs

**A-box weight:** 3.1 kg

**F-box dimensions:** 22.0 x 12.7 x 18.3

(HxWxD) cm

## Technical specifications

**Insulation:** Class II (double isolation)

**Type of lamps:** R95, 100 W

**Lifetime of lamp:** 750 sessions of 10 minutes

## Safety

**IEC certified:** Complies to IEC 60335

## Logistic data

**Country of origin:** Germany

## Medical appliance

**Medical Device Directive:** 2007/47/EC, MDD

93/42/EEC

## Logistic data

**CTV code:** 886151101000

