

PHILIPS

InfraCare infrared
lamp

650 W

Half body treatment

124 cm high

 Product Picture

HP3643



Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare provides effective pain relief for muscles and joints. Its 650 W infrared lamp is very comfortable and the infrared warmth penetrates deep into the skin, stimulates blood circulation and provides relief for half-body areas of 60x40 cm.

Easy to use

- Adjustable angle with rotation
- Easily reaches neck and shoulders

Effective pain relief

- Deep penetrating Infrared warmth

Half-body treatment

- 650 Watt infrared halogen lamp

Stimulated blood circulation

- Stimulates blood circulation

Highlights

650 Watt infrared halogen lamp



The Philips InfraCare HP 3643 has been optimised for effective relief of pains in areas covering 60 x 40 cm, such as the entire back, both shoulders and neck or thigh. The InfraCare features innovative infrared halogen lamp technology. The special optics, filter and the powerful halogen lamp of 650 Watt have been optimised, enabling half-body treatment (60 x 40 cm), and much more comfortable warmth, as the infrared warmth is more equally divided over the treatment area.

Adjustable angle with rotation

The InfraCare HP 3643 is easy to use for half-body treatment and for different situations (e.g. lying on a bed or sofa) due to its maximum flexibility. In addition to the height adjustment, the lamp house can easily be positioned for optimum treatment of the right body part by moving it into a horizontal or vertical position (left/right 90 degrees) and by tilting it up or down (45 degrees).

Easily reaches neck and shoulders

The InfraCare HP3643 can easily reach neck shoulders to give optimal warmth treatment of sore muscles. The height of the lamp house can be extended from 60 to 120 cm, so it offers a convenient position whether you are lying on a bed or sitting on a chair.

Stimulates blood circulation

Infrared light and warmth penetrates deep into your skin, stimulating blood circulation and dilating blood vessels. Furthermore, infrared light speeds up the metabolic processes and the removal of waste substances from the body.

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

Specifications

Half-body treatment

Treatment area: 60 x 40 cm

Infrared halogen lamp: 650 W

Effective pain relief

Deep penetrating warmth

Technical specifications

Cord length: 300 m

Power: 650 W

Frequency: 50–60 Hz

Insulation: Class II (double insulation)

Voltage: 220–230 V, Taiwan 110 V

Lifetime of lamps: 500 hour(s)

Weight and dimensions

Product dimensions: 29 x 124 x 29 cm

Product weight: 7.2 kg

F-box weight: 9.1 kg

Qnt. on Euro pallet: 24 pcs

F-box dimensions (= A-box): 34.5 x 60 x 39 cm

F-box dimensions: 34.5 x 60 x 39 cm

F-box weight (= A-box): 9.1 kg

Logistic data

CTV code: 884364301000

Country of origin: Germany

Safety

15 minutes auto shut-off

IEC certified: Complies to IEC 60601 and 60335

Easy to use

2 handgrips: for easy transport and positioning

Adjustability: 0–40 degrees (backwards)

On/off switch

Extendable height: 62–127 cm

Easily positionable lamphouse (left/right): 90

Easily positionable lamphouse (up/down): 45

Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

