



Philips Avent  
Single Electronic Breast  
Pump



SCF292/01

# Designed for comfort

Inspired by nature

Our unique pump is designed with an electronic memory that puts you in control by learning and continuing your personal pumping style.

**Benefits and features:**

- Learns and continues your pumping rhythm

**Helps you breastfeed longer**

- More milk, naturally

**Stimulate and express**

- Soft massage cushions trigger natural let-down, like baby

**Use three ways**

- By mains power, battery power and manual

# Specifications

## Country of origin

- England

## Development stages

- Stage: 0–6 months

## Power

- Voltage: 100-240 V

## What is included

- Breast milk container (125 ml/4 oz): 1 pcs

- Bottle stand/funnel cover: 2 pcs
- Extra soft Newborn Flow teat: 1 pcs
- Spare parts: 1 pcs
- Breast milk container (260 ml/9 oz): 1 pcs

## What is included

- Battery pack: 1 pcs
- Electronic Breast Pump: 1
- Mains power cord: 1
- Manual pump parts for when power is not available: 1

# Highlights

## Soft massage cushions



The patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and are intended to help stimulate a fast natural let down.

## Mains, battery and manual



In addition to mains power, this single electronic pump can also be operated using battery power or it can convert to a manual pump, for total flexibility.

## Learns and continues your pumping rhythm



Electronic memory learns and, at the touch of a button, continues your personal pumping rhythm.

## More milk, naturally



Expression allows breast milk to be available, even when you can't be there



Issue date 2021-02-12

Version: 9.3.1

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)