

**PHILIPS**

Infrared lamp

HP3614



# Highlights

## Philips — the light expert

Philips is one of the world's leading companies in light, with over 100 years of light experience. Philips has already sold more than 15 million Infrared lamps worldwide.

# Specifications

## Easy to use

**Adjustability:** 0–40 degrees (backwards)

## Weight and dimensions

**Product weight:** 0.5 kg

**Product dimensions:** 20.0 x 12.0 x 16.0 (H x W x D) cm

## Technical specifications

**Cord length:** 1.6 m

**Power:** 100 W

**Voltage:** 220/230 V

**Frequency:** 50 Hz

## Weight and dimensions

**A-box dimensions:** 22.4 x 26.2 x 58.8 (H x W x D) cm

**F-box weight:** 1 kg

**Qty. on Euro pallet:** 288 pcs

**No. of F-boxes in A-box:** 6 pcs

**A-box weight:** 3.1 kg

**F-box dimensions:** 21.1 x 12.6 x 19.2 (H x W x D) cm

## Technical specifications

**Insulation:** Class II (double isolation)

**Type of lamps:** R95, 100 W

**Lifetime of lamp:** 750 sessions of 10 minutes

## Safety

**IEC certified:** In line with IEC 60335

## Logistic data

**Country of origin:** Germany

## Logistic data

**CTV code:** 884361401000

## Easy storage

**Cord storage**

