

Philips Avent  
Digital baby thermometer  
set



SCH540/01

## Gentle and accurate

Combining the benefits of a digital thermometer and a soother thermometer.

### **Convenient for you and your baby**

- Includes Avent Freeflow Soother

### **Fast and reliable digital temperature readings**

- Clear measurement displays and battery status indicators
- Battery low indicator

### **Comfort for your child**

- Soft, flexible tip digital thermometer
- Choice to adapt temperature measurement method to suit baby
- Specially designed soother thermometer

**PHILIPS**  
**AVENT**

# Specifications

## Country of origin

- Made in China: Thermometers, Yes

## Power

- Operating time on battery: >100 hour(s)
- Battery type: LR41
- Number of batteries: 1 per thermometer

## Accessories

- User manual

## Technical specifications

- Operating temperature range: 90–109.9 °F
- Measurement time: 10 secs ±3 secs in water (bath), 25 secs ±10 secs oral and rectal, 60 secs ±10 secs



- underarm
- Accuracy: 95–107.6 +/-0.2 °F

## Convenience

- Automatic shut-off

## Country of origin

- Made in England: Soothers, Yes

## Convenience

- Flexible tip
- Suitable for newborns
- Clear measurement display
- Battery status indicator

# Highlights

## Choice of measurement method

The Philips Avent thermometer accurately and conveniently measures your baby's body temperature orally, axially and rectally.

## Soft, flexible tip digital thermometer

Soft, flexible tip digital thermometer


## Specially designed soother thermometer

Enables easy transition from soother to temperature measurement.

## Clear measurement displays

Clear measurement displays and battery status indicators

## Battery low indicator



Battery low indicator

## Avent Freeflow Soother

Lets your baby get used to the shape.



Issue date 2021-02-12

Version: 5.1.1

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)