



Philips
Infrared muscular pain
relief

200 W
Focused treatment area



HP3621/02

Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 200 Watt infrared warmth is very comfortable and penetrates deep into skin, improving blood circulation thereby giving effective relief for focused areas of 20 x 30 cm.

PHILIPS

Specifications

Easy to use

- Adjustability: 0-40 degrees (backwards)
- On/off switch
- Soft-touch handgrip: for easy transport and positioning

Weight and dimensions

- Product dimensions: 20 x 30.9 x 17.2 (WxHxD) cm
- Product weight: 1.4 kg

Technical specifications

- Power: 200 W
- Cord length: 2 m
- Voltage: 110 or 220 or 230 V
- Frequency: 50/60 Hz

Weight and dimensions

- A-box dimensions: 29 x 40.4 x 74.5 (WxHxD) cm
- F-box weight: 1.8 kg
- No. of F-boxes in A-box: 3 pcs
- Qnt. on Euro pallet: 60 pcs
- A-box weight: 5.4 kg
- F-box dimensions: 28 x 38.6 x 24.5 (WxHxD) cm

Focused treatment

- Infrared halogen lamp: 200 W
- Treatment area: 20 x 30 cm

Technical specifications

- Insulation: Class II (double isolation)
- Lifetime of lamps: 500 hour(s)

Safety

- IEC certified: In line with IEC 60601 and 60335

Logistic data

- Country of origin: Germany

Logistic data

- CTV code: 884362102000

Effective pain relief

- Deep penetrating warmth
- Even warmth distribution: Now up to 3 x more comfortable warmth

Medical appliance

- Proven health benefits: supported by Medical Device Directive certificate

