

PHILIPS

Infrared muscular
pain relief

650 W

Half body treatment

63.5 cm high

HP3641



Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 650 W infrared warmth is very comfortable and penetrates deep into the skin, improving blood circulation and thereby giving effective relief for half-body areas of 60 x 40 cm.

Easy to use

- Adjustable angle with rotation

Proven health benefits


- Supported by the Medical Device Directive (MDD)

Treatment control

- Digital timer with automatic shut-off

Highlights

Adjustable angle with rotation



The InfraCare HP 3643 is easy to use for half-body treatment and for different situations (e.g. lying on a bed or sofa) due to its maximum flexibility. In addition to the height adjustment, the lamp house can easily be positioned for optimum treatment of the right body part by moving it into a horizontal or vertical position (left/right 90 degrees) and by tilting it up or down (45 degrees).

Digital timer with automatic shut-off



Treatment time can easily be controlled by setting the advised or preferred number of minutes. The digital timer automatically counts down and automatically shuts off at the end of the session.

Supported by the Medical Device Directive (MDD)



This device has proven health benefits, as certified by the Medical Device Directive (MDD).

Specifications

Easy to use

2 handgrips: for easy transport and positioning

On/off switch

Digital user interface: for feedback on time and warmth intensity

Extendable height: 45-80 cm

Easily positionable lamphouse (left/right): 90

Easily positionable lamphouse (up/down): 45

Weight and dimensions

Product dimensions: 29.0 x 63.5 x 29.0 (WxHxD) cm

Product weight: 5.5 kg

Technical specifications

Power: 650 W

Cord length: 3 m

Voltage: 220 V

Frequency: 50 and 60 Hz

Weight and dimensions

Qty. on Euro pallet: 18 pcs

F-box dimensions (= A-box): 33.9 x 80.0 x 36.5 (WxHxD) cm

F-box weight (= A-box): 6.5 kg

Technical specifications

Insulation: Class II (double isolation)

Lifetime of lamps: 500 hour(s)

Half-body treatment

Infrared halogen lamp: 650 W

Treatment area: 60 x 40 cm

Safety

IEC certified: In line with IEC 60601 and 60335

Logistic data

Country of origin: Germany

Logistic data

CTV code: 884364101000

Effective pain relief

Deep penetrating warmth

Even warmth distribution: Now up to 3 x more comfortable warmth

Medical appliance

Proven health benefits: supported by Medical Device Directive certificate

Treatment control

Digital timer: with auto shut-off

Dimmer with 3 warmth settings

Pause function

