



Philips  
Infrared lamp



HP3614

## Relieves muscular pains and common colds

### Deep, focused warmth

Philips infrared lamp gives pain relief for muscles and joints. Its 100 Watt infrared warmth is comfortable and penetrates deep into skin, improving blood circulation thereby giving relief for focused areas of 20x30 cm.

#### **Philips - the light expert**

- More than 15 million infrared lamps sold

# PHILIPS

# Specifications

## Easy to use

- Adjustability: 0-40 degrees (backwards)

## Weight and dimensions

- Product weight: 0.5 kg
- Product dimensions: 20.0x12.0x16.0 (HxWxD) cm

## Technical specifications

- Cord length: 1.6 m
- Power: 100 W
- Voltage: 220/230 V
- Frequency: 50 Hz

## Weight and dimensions

- A-box dimensions: 22.4x26.2x58.8 (HxWxD) cm
- F-box weight: 1 kg
- A-box weight: 3.1 kg
- Qnt. on Euro pallet: 288 pcs

- No. of F-boxes in A-box: 6 pcs
- F-box dimensions: 21.1x12.6x19.2 (HxWxD) cm

## Technical specifications

- Insulation: Class II (double isolation)
- Type of lamps: R95, 100W
- Lifetime of lamp: 750 sessions of 10 minutes

## Safety

- IEC certified: In line with IEC 60335

## Logistic data

- Country of origin: Germany

## Logistic data

- CTV code: 884361401000

## Easy storage

- Cord storage

# Highlights

## Philips - the light expert

Philips is one of the global leading companies in light, with over 100 years of light experience. Globally Philips already sold more than 15 million of Infrared lamps.



Issue date 2021-02-15

Version: 3.1.1

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)