

Philips AVENT  
Disposable breast pads

2 day pads

SCF254/02



## Comfortable and discreet

### AVENT breast pads

Unique Philips AVENT breast pads specially designed to help you stay dry and comfortable during the day.

#### **Feel confident and relaxed**

- Invisible under clothes
- Anti slip with adhesive tapes
- Developed with breastfeeding expert
- Hygienic
- Silky soft feel
- Soft with discreet, contoured shape

#### **No leakage of breast milk**

- Unique four layered construction for maximum dryness

**PHILIPS**  
**AVENT**

Disposable breast pads  
2 day pads

# Specifications

## Country of origin

- China

## Material

- Breast pads: Dermatologically tested

## Design

- Discreet contoured shape

## What is included

- Disposable day breast pads: 2 pcs

## Dimensions & weight

- Dimensions: 100x25 mm

## Functions

- Ultra dry: Multi layers, One-way top

## Maximum comfort

- Anti-slip: Adhesive tape

## Silky soft feel

- Silky soft topsheet: Breathable natural material, Dermatologically tested

SCF254/02

# Highlights

## Invisible under clothes



Ultra thin with a contoured shape. Adhesive strip to keep the pads in place.

## Anti slip adhesive tape

Anti slip due to adhesive tapes to keep the breast pad in place.

## Developed with expert

Developed in collaboration with a midwife and breastfeeding advisor who has been supporting mums with breastfeeding for 20 years.

## Hygienic

Individually wrapped for extra hygiene.

## Silky soft feel



Silky soft topsheet and breathable, natural materials. Dermatologically tested.

## Soft with discreet, contoured shape

The Philips Avent breast pads have a contoured shape that allows you to discretely wearing under your clothing.

## Unique, layered construction for maximum dryness

1. Ultra soft top layer with nipple indent- keeps breast dry at all times. 2. Ultra absorbent layer, cushioned for extra comfort. 3. Ultra absorbent core- draws in moisture to guard against embarrassment. 4. Ultra breathable outer layer- helps avoid sore nipples. Anti-slip with adhesive tapes- keeps pad in place.



Issue date 2021-02-15

Version: 7.2.1

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)