



Philips Avent  
Manual breast pump

with VIA Storage System

SCF300/12

## Express more. Store easily.

Express, store and feed in easy steps

The unique Philips Avent Manual Breast Pump SCF300/12 offers you total control of the suction level and pumping rhythm. It assures you enhanced efficiency, being clinically proven to express more milk than a hospital electric pump.

**Benefits and features:**

- Gently stimulates fast milk flow

**Clinically proven effective**

- Patented Let-down Massage Cushion

**Fits all Philips Avent breast pumps and teats**

- Use adaptor to express and feed

**For safe storage and transport**

- Leak-proof twist-on system

**PHILIPS**  
**AVENT**

# Specifications

## Country of origin

- England

## Development stages

- Stage: 0 - 6 months

## Material

- BPA free\*

## What is included

- Milk storage cup lid: 5 pcs
- Milk storage cup adaptor: 2 pcs
- Newborn nipple travel pack: 2 pcs
- Extra soft Newborn Flow nipple: 2 pcs
- Spare parts: 1 pcs
- Pre-sterilized VIA Cup (180 ml/ 6 oz): 5 pcs
- Manual breast pump: 1 pcs

SCF300/12

# Highlights

## Gently stimulates fast milk flow

As you pump the patented Philips Avent Let-down Massage Cushion flexes in and out, gently massaging the area around the nipple to help stimulate natural milk flow.

## Patented Let-down Massage Cushion



Five, soft petal massagers of the Philips Avent breast pump gently flex in and out as you pump and work with the vacuum to imitate your baby's suckling.

## Leak-proof twist-on system



For safe storage and transport

## Adaptor to express and feed



Use the adaptor to express and feed with the same cup.



Issue date 2021-02-15

Version: 6.0.1

© 2021 Koninklijke Philips N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips N.V.  
or their respective owners.

[www.philips.com](http://www.philips.com)

\* 0% BPA, following EU regulation 10/2011