

PHILIPS

 Product Picture

AVENT

Single electronic
breast pump

SCF292/01

Designed for Comfort

Inspired by nature

Our unique pump is designed with an electronic memory that puts you in control by learning and continuing your personal pumping style.

Benefits and features:

- Learns and continues your pumping rhythm

Helps you breastfeed longer

- More milk, naturally

Stimulate and Express

- Soft massage cushions trigger natural let-down, like baby

Use three ways

- By mains power, battery power and manual

Highlights

Soft massage cushions



Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

Mains, battery and manual

In addition to mains power, this single electronic pump can also operate using battery power and as a manual pump, for total flexibility.

Learns and continues your pumping rhythm

Electronic memory learns and at the touch of a button, continues your personal pumping rhythm.

More milk, naturally



Expression allows breast milk to be available, even when you can't be there

Specifications

Country of origin

England

Development stages

Stage: 0 - 6 months

Power

Voltage: 100 - 240 V

What is included

- Breast milk container (125 ml/ 4 oz): 1 pcs
- Bottle stand/funnel cover: 2 pcs
- Extra soft Newborn Flow nipple: 1 pcs
- Spare parts: 1 pcs
- Breast milk container (260 ml/ 9 oz): 1 pcs

What is included

- Electronic Breast Pump: 1
- Battery pack: 1 pcs
- Mains power lead: 1
- Manual pump parts for when power is not available: 1

