

Philips Avent
Twin electronic breast
pump



SCF304/02

Designed for comfort

Inspired by nature

Our unique BPA-Free pump is designed with an electronic memory that puts you in control by learning and continuing your personal pumping style.

Benefits and features:

- Ultra comfortable and effective
- Control at your fingertips
- Infinitely variable speed, vacuum and rhythm – no pre-sets
- Quiet, effortless and portable
- Pumps more milk faster

Stimulate and Express

- Soft massage cushions trigger natural let-down, like baby
- Gentle vacuum mimics baby's suckling for steady milk flow

PHILIPS
AVENT

Highlights

Gentle-draw vacuum

The gentle-draw vacuum mimics baby's suckling to create a steady milk flow that requires less pumping.

Soft massage cushions



Patented massage cushions of the Philips Avent breast pump flex in and out, replicating baby's suckling action and is intended to help stimulate a fast natural let down.

Infinitely variable control

Enjoy complete control over speed, vacuum and pattern – you can even set the pump to match your baby's own suckling action.

Intuitive control handle



Everything that controls the pump is at your fingertips and close to the breast, which means you can make adjustments instantly and intuitively- no need to reach elsewhere.

Patented Let-down Massage Cushion



Five, soft petal massagers gently flex in and out working in tandem with the reliable vacuum. This is intended to help stimulate a fast, natural milk flow.

Pumps more milk faster



Double pumping has been shown to improve milk production. The twin electronic breast pump is the only breast pump designed to maintain the same level of vacuum and control for both pumps, which helps to ensure effective use.

Quiet, effortless and portable



Quiet, effortless and portable

Specifications

Country of origin

- England

Development stages

- Stage: 0 - 6 months

Power

- Voltage: 100 - 240 V

What is included

- Bottle stand/funnel cover: 4 pcs
- Breast milk container (125 ml/ 4 oz): 2 pcs
- Extra soft newborn flow nipple: 2 pcs
- Sealing disc for milk storage: 2 pcs
- Newborn nipple travel pack: 2 pcs
- Cool pack: 4 pcs
- Spare parts: 1 pcs

What is included

- Twin Electronic Breast pump: 1 pcs
- Microfiber travel bag: 1 pcs
- Thinsulate™ carrier: 2 pcs



Issue date 2021-02-15

Version: 9.0.1

© 2021 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com